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---

### Legend

**(Representation of individual parts: see Fig. A in enclosed leaflet.)**

- 1 Headgear YARAstrap
- 2 Forehead cushion
- 3 Fine adjustment component of forehead support
- 4 Coarse adjustment component of forehead support
- 5 Forehead support
- 6 Headgear clip
- 7 Pressure measurement ports
- 8 Port cap
- 9 Anti-asphyxia valve (YARA only)
- 10 Rotating sleeve\*
- 11 Elbow\*
- 12 Retaining ring\*
- 13 Mask frame
- 14 Mask cushion

\*YARA: transparent; YARA NV: blue

### Intended use

The YARA and YARA NV full face masks are accessories. They act as a sealing element between the patient and the therapy unit, and are used for the treatment of sleep apnea and for the non-invasive ventilation of patients with ventilatory insufficiency.

The YARA NV mask must only be used in combination with ventilation devices that have an active exhalation valve.

### Before use by the patient

Patients must be instructed in the use of the mask by qualified personnel before therapy begins.

## 2 Safety instructions

---

### Safety regulations

Following the requirements of Directive 93/42/EEC please observe the following points:

- Please read these operating instructions through carefully. They are a part of the product and should be available at all times.
- Never use YARA if the anti-asphyxia valve is damaged or missing, as this may lead to increased CO<sub>2</sub> reinhalation if the therapy unit fails.
- Never use YARA or YARA NV for patients who are not able to remove it independently or have a tendency to vomit. There is a risk of asphyxia.
- Only use the YARA or YARA NV at therapy pressures between 4 - 30 hPa or 4 - 40 hPa (models YARA 40 hPa / YARA NV 40 hPa). Otherwise it can lead to CO<sub>2</sub> reinhalation or to leakages.
- Malfunctions and a lack of biocompatibility may result if third-party articles are used. Please bear in mind that in these cases any warranty entitlement and liability shall lapse where the accessories recommended in the instructions for use or original spare parts are not utilized.
- Patients with limited spontaneous respiration must be monitored at all times. The appropriate negative pressure/leak alarm on the ventilation device must be activated.
- Patients with obstructive and restrictive lung diseases should be examined on an individual basis during use of the YARA or YARA NV (e.g. with regard to the reinhalation of CO<sub>2</sub>).
- Never close the opening of the anti-asphyxia valve or the exhaled air system. Otherwise there is a risk of reinhaling CO<sub>2</sub> and of asphyxia.

- Do not place the mask over your face for a prolonged period unless the therapy unit is switched on. Otherwise there is the risk of reinhaling CO<sub>2</sub>.
- Do not use a breathing hose that is electrically conductive.
- Follow the instructions in section 4 "Cleaning and disinfecting instructions" to avoid infection, bacterial contamination and poor functionality of the mask.

### Contraindications:

Under certain circumstances the YARA or YARA NV should not be used or only be used with particular care. In each case the decision whether to use these masks for therapy lies with the doctor in charge.

- Functional disturbance of the lower oesophageal sphincter (lower sphincter of the oesophagus)
- Excessive regurgitation (reflux of stomach contents into the oesophagus)
- Cough reflex disturbance
- Hiatus hernia (diaphragmatic hernia)
- Anesthetic or sedation
- Open facial wounds which tend towards inflammation
- Haemodynamic or cardiorespiratory instability
- Unconsciousness
- Claustrophobia, anxiety
- Facial deformation or deformation of the end of the nasopharynx
- Beard or other obstructions to adjusting and sealing the mask
- Facial injuries
- Barotraumas
- Necessity for ventilation or artificial respiration 24 hours a day
- Fresh facial, oesophageal or stomach operations
- Taking medication that can cause vomiting
- Necessity for immediate intubation

### Side effects

The following side effects may occur:

Nasal congestion, dry nose, dry mouth in the morning, sinus pressure, chest pains, breathlessness, flatulence, headaches, irritation of the conjunctiva, irritating noises during exhalation, pressure sores on the face, reddening of the skin.

## 3 Operation

---

**Fig. A-1: See enclosed leaflet.**

### 3.1 Adjust forehead support (Fig. C)

#### Fine adjustment component

Squeeze the upper side tabs of the coarse adjustment component inwards. Position the fine adjustment component in one of the three slots.

To adjust the height of the forehead cushion or the fine adjustment component (approximately 1 cm), turn the fine adjustment component through 180°. Then position this as described above.

#### Coarse adjustment component

Squeeze the lower side tabs of the coarse adjustment component inwards. Place them in one of the three slots in the mask frame.

### 3.2 Put on the mask (Fig. B)

#### Caution!

To minimize the risk of vomiting, the patient should not eat or drink anything for 3 hours before using the mask.

1. Unhook one of the lower headgear clips from the mask.
2. Place the headgear over your head and hook the headgear clip back onto the strap.
3. Position the mask on your face as shown in Figure B.
4. Adjust the upper and lower headgear straps so that the mask sits firmly, but not too tightly, on your face. If necessary, re-adjust the forehead support.

### 3.3 Exhaled air system (Fig. H)

#### YARA

Yara (all mask components are transparent) has an integrated exhaled air system. Do **not** connect an external exhaled air system. Otherwise the therapy pressure will not be reached.

#### YARA NV

YARA NV (retaining ring, elbow and rotating sleeve colored blue) has **no** exhaled air system.

YARA NV must **only** be used in combination with ventilation devices that have an **active exhalation valve**. Follow the instructions for use of your ventilation device.

If a therapy unit is used without an active exhalation valve, there is a risk of asphyxia should the unit fail.

Do not connect **an external passive** exhaled air system (for example Silentflow 2). Otherwise smooth functioning of the therapy unit cannot be guaranteed.

### 3.4 Anti-asphyxia valve (Fig. G)

**Only YARA** is equipped with an anti-asphyxia valve. If the therapy unit should fail, this valve allows the patient to breathe the ambient air. The risk of reinhalation of CO<sub>2</sub> is significantly reduced and thus also the risk of asphyxia.

#### Caution!

The anti-asphyxia valve is a very sensitive component. Always handle it with care. Always check that the anti-asphyxia valve is working before each use as described in section 5 "Functional check".

### 3.5 Connecting the mask to the therapy unit (Fig. H)

Push the hose system of your therapy unit onto the rotating sleeve. Ensure that the connection to the hose is secure and airtight. The connection must not be loosened by the therapy pressure.

If you want to detach the mask temporarily from the hose, pull the rotating sleeve off the elbow.

#### Caution!

If there is a high risk of transferring bacteria between the patient and the unit, a bacterial filter should be used. Also pay attention to the instructions for use of the therapy unit.

### 3.6 Ports for pressure measurement (Fig. E)

YARA and YARA NV are fitted with two ports that can be used to measure pressure or to feed in oxygen.

#### Caution!

- Oxygen therapy can have side effects if incorrectly administered or used in the wrong dosage. Only carry out this sort of treatment if prescribed by a doctor.
- When carrying out oxygen therapy, strictly follow the instructions for use of your O<sub>2</sub> system and your therapy unit.

When the pressure measurement ports are not in use, always close them with the caps provided. Otherwise air will escape from the ports and the therapy pressure will not be reached.

### 3.7 Compatible devices and accessories

#### Therapy units

YARA and YARA NV masks can be combined with all units designed for carrying out the therapies described under "Intended Use".

Note that if unsuitable device combinations are used, the actual pressure in the mask may not correspond to the therapy pressure you have been prescribed.

This may also happen even if the correct pressure is indicated on the therapy unit. Have your device combination adjusted by a doctor or distributor so that the pressure in the mask corresponds to your therapy pressure.

#### Humidifier

YARA and YARA NV can be used with both cold and warm-air humidifiers.

Never set the humidification level so high that condensation forms in the air hose. Follow the instructions for use of your humidifier and therapy unit.

#### Quick release system (optional) (Fig. I)

To open the headgear in an emergency (e.g. vomiting) a quick release system can be supplied as an accessory.

The quick release system consists of a special clip with ripcord, which is fitted instead of one of the lower headgear clips.

In an emergency, simply pull the ripcord. The special clip will detach itself from the mask. The headgear is now open and the mask can be removed.

### 3.8 After usage

1. Undo one of the lower headgear clips from the mask and remove the mask.
2. Undo the other three headgear clips in order to remove the headgear from the mask.
3. Proceed as described in "4 Cleaning and disinfecting instructions".

### 3.9 Disassembling/assembling the mask

#### Disassembling the mask (Fig. G - C)

1. Detach the mask cushion from the mask frame.
2. Remove the cap from the pressure measurement ports.
3. Turn the retaining ring anticlockwise to the 11 o'clock position. Now pull the retaining ring and the elbow away from the mask frame.
4. Press the elbow out of the retaining ring and pull the rotating sleeve off the elbow.
5. **Only YARA:** Remove the anti-asphyxia valve carefully off the elbow. Ensure that the membrane on the interior of the valve is not damaged.
6. Detach the coarse adjustment component from the mask frame.
7. Detach the fine adjustment component from the coarse adjustment component. Then detach the forehead cushion from the fine adjustment component.

#### Assembling the mask (Fig. C - G)

1. Attach the forehead cushion to the fine adjustment component. Then attach the fine adjustment component to the coarse adjustment component.
2. Attach the coarse adjustment component to the mask frame.
3. Push the rotating sleeve onto the elbow until you hear it click into place.
4. **Only YARA:** Push the anti-asphyxia valve carefully onto the elbow until you hear it click into place on both sides of the elbow. Take care that the membrane on the interior of the valve is not damaged or jammed.
5. Press the rotating sleeve onto the elbow. Then press the elbow into the retaining ring until you hear it click into place.

6. Push the retaining ring (in the 11 o'clock position) and elbow onto the mask frame. Turn the retaining ring to the 12 o'clock position.
7. Push the cap onto the pressure measurement ports.
8. Attach the mask cushion to the mask frame.
9. Fasten the headgear to the mask with the clips.
10. Check that everything is working properly. See "Functional check" on page 17.

## 4 Cleaning and disinfecting instructions

### 4.1 Intervals

#### Daily

Clean the mask in warm water to remove all residues. To do this the mask can be disassembled as described in section 3.9 "Disassembling/assembling the mask". Use a mild cleaning agent (e.g. dishwasher detergent). After cleaning, rinse all parts carefully with clean water.

#### Caution!

For patients with a weakened immune system or a special medical history, **daily** disinfecting or sterilization of the mask components may be necessary. In this case, please consult your doctor.

#### Weekly

Disassemble the mask as described in section 3.9 "Disassembling/assembling the mask" and clean and sterilize the parts of the mask. See "Cleaning, disinfection, sterilization" on page 16.

## 4.2 Cleaning, disinfection, sterilization

### Permissible procedures

Procedure	Dishwasher	Washing machine	Disinfection*	Sterilization (alternative to disinfection)
Mask cushion	•		•	•
Forehead cushion	•		•	
Port cap	•		•	•
Mask frame	•		•	
Coarse adjustment component	•		•	
Fine adjustment component	•		•	
Elbow	•		•	
Anti-asphyxia valve			•	
Rotating sleeve	•		•	
Retaining ring	•		•	
Headgear clip	•	•	•	
Headgear		•		
Quick release system		•		

#### Caution!

- Always clean the components of the mask thoroughly before disinfecting or sterilizing them as indicated in the above table
- The anti-asphyxia valve is a very sensitive component. Always handle it with care to avoid damage or deformation.

#### Description of procedures

Dishwasher	Rinse cycle at 65°C. Then carefully rinse all parts with clean water.
Washing machine	Rinse cycle at 30°C.
Disinfection*	Clean mechanically, e.g. with a brush, in dilute solution. After disinfection rinse all parts thoroughly with distilled water and dry thoroughly.
Sterilization (alternative to disinfection)	Steam sterilization in devices to comply with EN 285. Temperature 134 °C, minimum hold time 3 minutes.

\*We recommend the disinfectant GIGASEPT FF. It has been proven that this disinfectant is compatible. Follow the manufacturer's instructions without fail.

### Note

We recommend that you wear gloves (e.g. disposable gloves) while cleaning and disinfecting your mask.

## 4.3 After cleaning and disinfecting

Dry all parts thoroughly after cleaning and disinfecting. Check the components for residual contamination or damage. If necessary, repeat the cleaning and disinfecting. Replace damaged parts. Reassemble the mask.

## 4.4 Multiple patient use

If the mask is used by a different patient, it must be cleaned and disinfected first.

Disassemble the mask. Disinfect or sterilize all parts as described in section "4.2 Cleaning, disinfection, sterilization".

Use a new headgear and where necessary a new quick release system.

4. Switch on the therapy unit. The membrane should now close the aperture of the anti-asphyxia valve. Only then can ambient air flow from the therapy unit into the mask.

5. Do not use the mask if the valve or the membrane is not working. Replace the anti-asphyxia valve. See "Scope of supply" on page 18.

### YARA and YARA NV

Carry out a visual check after each cleaning and disinfecting procedure.

If the mask cushion is damaged and is no longer airtight, replace it with a new one.

Check the functionality and airtightness of the ball-and-socket joint on the elbow. If it no longer moves easily or is no longer airtight, replace the elbow and retaining ring.

Parts that are no longer usable can be disposed of as normal household waste.

## 5 Functional check

---

Before use check that everything is working properly as described below.

### YARA only

Always check that the anti-asphyxia valve is working before each use as described below:

1. Ensure that the valve and membrane are dry, undamaged and free of secretions. Do not use the mask if the valve or membrane is damaged (e.g. through tears and deformation). In this case replace the anti-asphyxia valve. See "Scope of supply" on page 18.
2. Ensure that the valve membrane is not closing the valve aperture. Only then can ambient air flow into the mask.
3. Attach the breathing hose (delivered with the therapy unit) to the elbow of the mask and to the therapy unit.

## 6 Troubleshooting

Mal-functions	Cause	Solution
Painful pressure on the face.	Mask is too tight.	Loosen headgear slightly. Adjust position of forehead support.
Draught in the eyes.	Mask is too loose.	Tighten headgear slightly. Adjust position of forehead support.
	Mask does not fit.	Use mask of different size. Available from your distributor.
	High therapy pressure.	Use hard mask cushion (for pressures up to 40 hPa). Available from your distributor.
Therapy pressure is not reached.	Mask not adjusted correctly.	Adjust headgear. Adjust position of forehead support.
	Mask cushion is damaged.	Replace mask cushion.
	Hose system is not airtight.	Check push-fit connector and fit of hose sleeves.
	Air is escaping from pressure measurement ports.	Cover pressure measurement ports with cap (included in package).
	Anti-asphyxia valve: – incorrectly fitted – defective	Anti-asphyxia valve:  fit correctly  replace

## 7 Scope of supply

### Note

All masks are supplied fully assembled (including headgear YARAstrap).

### YARA models

Article	Article no.
YARA, size S	WM 26410
YARA, size M	WM 26420
YARA, size L	WM 26430
YARA 40 hPa, size S	WM 26411
YARA 40 hPa, size M	WM 26421
YARA 40 hPa, size L	WM 26431

### YARA NV models

Article	Article no.
YARA NV 40 hPa, size S	WM 26461
YARA NV 40 hPa, size M	WM 26471
YARA NV 40 hPa, size L	WM 26481

### Replacement parts

Article	Article no.
Ball-and-socket joint set: Retaining ring, elbow, rotating sleeve, anti-asphyxia valve	WM 26569
Ball-and-socket joint set NV (blue): Retaining ring NV, elbow NV, rotating sleeve NV	WM 26207
Forehead support: Forehead cushion, coarse adjustment component forehead support, fine adjustment component forehead support	WM 26201
Forehead cushion	WM 26200

Article	Article no.
Mask cushion 30 hPa, size S	WM 26510
Mask cushion 30 hPa, size M	WM 26520
Mask cushion 30 hPa, size L	WM 26530
Mask cushion 40 hPa, size S	WM 26511
Mask cushion 40 hPa, size M	WM 26521
Mask cushion 40 hPa, size L	WM 26531
Headgear, blue: YARAstrap, inc. headgear clips	WM 26290
Headgear clips (4 pcs.)	WM 26205
Port cap	WM 26208
Instructions for use DE, GB, FR, NL, IT	WM 66117

### Accessories

Article	Article no.
Hair protector, blue: YARAcap, inc. headgear clips	WM 26293
Quick release system (ripcord)	WM 26560

## 8 Specifications

### Technical data

Parameters	Value
Device class according to Directive 93/42 ECC:	II a
Dimensions (WxHxD):	c. 107 x 160 x 106 mm
Weight:	c. 102 grams
Therapy pressure range: Mask cushion 30 hPa: Mask cushion 40 hPa:	4 to 30 hPa 4 to 40 hPa
Hose connector: cone complying with EN ISO 5356-1 YARA: YARA NV:	∅ 22 mm (male) ∅ 22 mm (female)

Parameters	Value
Pressure measurement ports:	∅ 4 mm
Width of headgear:	max. 22 mm
Temperature range Operation: Storage	+ 5° C to + 40° C - 20° C to + 70° C
Flow resistance YARA: at 50 l/min at 100 l/min YARA NV: at 50 l/min at 100 l/min Anti-asphyxia valve (YARA) at 50 l/min at 100 l/min	0.02 hPa 0.26 hPa 0.09 hPa 0.37 hPa 0.5 hPa 2.0 hPa
Sound pressure level at 10 hPa	27.5 dB(A)
Switching pressure anti- asphyxia valve (YARA) Open: Close:	≤ 2 hPa ≤ 1 hPa
Working life	6 to 12 months*

\* depending on frequency of use and cleaning

### Materials


Article	Material
Mask cushion, port cap	Silicone
Forehead cushion, retaining ring	TPE, Polypropylene
Mask frame, elbow	Polycarbonate
Headgear YARAstrap	Cotton, polyamide, elastane
Rotating sleeve, forehead support coarse adjustment component, forehead support fine adjustment component	Polypropylene
Headgear clip	POM

Article	Material
Valve cartridge (anti-asphyxia valve)	TPE, Polypropylene
Quick release system (ripcord)	POM, Polyester, Polyamide

All mask components are latex-free.

### Pressure/flow characteristic curve

The pressure/flow characteristic curve shows outlet flow as a function of therapy pressure.

(Subject to design changes)

## 9 Warranty

- Weinmann guarantees that the product is free from defects, if used as intended, for a period of 6 months from date of purchase.
- Warranty claims must be accompanied by proof of purchase showing merchant and date of purchase.
- We acknowledge no liability under warranty in the event of:
  - Non-compliance with the instructions for use
  - Incorrect operation
  - Improper use or treatment
  - Repairs to the device attempted by unauthorized persons
  - Force majeure
  - Transport damage due to inadequate packing of returned items
  - Deterioration in normal use and customary wear
  - Failure to use original replacement parts.

- Weinmann does not accept liability for damage consequential upon defects unless due to willful acts or gross negligence or in the event of harm to life or limb due to minor negligence.
- Weinmann reserves the right to decide whether to eliminate defects, to deliver a defect-free item or to reduce the purchase price by a reasonable amount.
- If a warranty claim is refused, we will not accept the transport cost to or from our address.
- This does not affect your statutory rights.

EEC Letter of Conformity on Medical Devices	
We, the manufacturer	Weinmann Geräte für Medizin GmbH+Co.KG Kronsassweg 40 · 22525 Hamburg – Germany
declare in sole responsibility the	Product Name: <b>Full Face Mask</b> Type/Model: <b>YARA size S</b> <b>YARA size M</b> <b>YARA size L</b>
being in conformity with the respective regulations of the following guideline:	Medical Device Directive 93/42/EWG
Classification:	IIa
standards used	in particular: DIN EN ISO 17510-2 (May 2003)
Mark:	The product has the mark by TUV Rheinland Product Safety GmbH 51105 Cologne – Germany
Hamburg, 07./03./2006	
 Dr. K.-A. Feldhahn President Chief Technology Officer	
<b>WEINMANN</b>	
<small>www.71071.07.06</small>	

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Il prodotto è provvisto di marchatura del TÜV Rheinland Product Safety GmbH
Geräte für Medizin GmbH+Co. KG