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1 Device description

Legend

(Description of individual parts: see enclosed leaflet)

- 1 Headgear JOYCE® *strap*
- 2 Forehead cushion
- 3 Fine adjustment component of forehead support
- 4 Coarse adjustment component of forehead support
- 5 Forehead support
- 6 Headgear clips
- 7 Pressure measurement ports
- 8 Port cap
- 9 Rotating sleeve (NV: matt finish)
- 10 Elbow (NV: matt finish)
- 11 Retaining ring (NV: matt finish)
- 12 Mask frame
- 13 Mask cushion

Intended use

The nasal mask is an accessory and acts as a seal between the patient and the therapy device. It is used for the treatment of sleep apnoea and for the non-invasive ventilation of patients with ventilatory insufficiency.

Before use by the patient

Patients must be instructed in the use of the mask by qualified personnel before therapy begins.

2 Safety instructions

Safety regulations

Following the requirements of Directive 93/42 ECC please observe the following points:

- Please read these instructions for use very carefully. They are a part of the product and should be available at all times.
- Use the mask only with therapy pressures between 4 and 30 hPa or between 4 and 40 hPa (other models JOYCE® 40 hPa / JOYCE® NV 40 hPa).
- Malfunctions and a lack of biocompatibility may result if third-party articles are used. Please bear in mind that in these cases any warranty entitlement and liability shall lapse where the accessories recommended in the instructions for use or original spare parts are not utilized.
- Patients with limited spontaneous respiration must be monitored at all times. The appropriate negative pressure/leak alarm on the ventilation device must be activated.
- Patients with obstructive and restrictive lung diseases should be examined on an individual basis during use of the JOYCE® (for example with regard to re-inhalation of CO₂).
- Do not close the openings in this mask except in situations described in these instructions for use. Otherwise there is a risk of re-inhaling CO₂ and of suffocation.
- Only put the nasal mask over your face for a prolonged period when the therapy device is switched on. Otherwise there is a risk of re-inhaling CO₂.
- Do not use a breathing hose that is electrically conductive.
- Follow the instructions in section 4 "Cleaning and disinfecting instructions" to avoid infection, bacterial contamination and poor functionality of the mask.

Side-effects

The following side-effects may occur:

nasal congestion, dry nose, dry mouth in the morning, sinus pressure, irritation of the conjunctiva, irritating noises during exhalation, pressure sores on the face, reddening of the skin.

3 Operation

Fig. A-E: see enclosed leaflet.

3.1 Adjust forehead support (Fig. A)

Fine adjustment

Squeeze the upper side tabs of the coarse adjustment component inwards. Fit the fine adjustment component into one of the three notches.

To change the height of the forehead cushion, i.e. the fine adjustment (by approx. 1 cm), turn the fine adjustment component through 180°. Then position it as described above.

Coarse adjustment

Squeeze the lower side tabs of the coarse adjustment component inwards. Place it in one of the slots in the mask frame.

3.2 Putting on the mask

1. Unhook one of the lower headgear clips from the mask.
2. Guide the headgear over your head and hook the headgear clips back onto the strap.
3. Adjust the upper and lower headgear straps so that the mask sits firmly, but not too tightly, on your face.
4. If necessary, re-adjust the forehead support.

3.3 Exhaled air system (Fig. E)

JOYCE® (Vented)

In the vented model (all components of the mask have a polished surface) JOYCE® has an integral exhaled air system. Do **not** connect an external exhaled air system. Otherwise the therapy pressure will not be reached.

JOYCE® NV (Non-Vented)

In the Non-Vented model (retaining ring, elbow and rotating sleeve with matt surfaces) JOYCE® does not have an exhaled air system. An external exhaled air system must be connected **without exception**. Otherwise

there is a risk of re-inhaling CO₂ and of suffocation.

Follow the instructions for use of your exhaled air system.

3.4 Connecting the mask to the therapy device (Fig. E)

Push the hose of your therapy device onto the rotating sleeve (Vented) or onto your external exhaled air system (Non-Vented). Ensure that the connection to the hose is secure and airtight. The connection must not be loosened by the therapy pressure.

If you want to detach the mask temporarily from the hose, pull the rotating sleeve off the elbow.

3.5 Ports for pressure measurement (Fig. C)

JOYCE® masks are fitted with two ports that can be used to measure pressure or to feed in oxygen.

Caution!

- Oxygen therapy can have side-effects if incorrectly administered or used in the wrong dosage. Only carry out this sort of treatment if prescribed by a physician.
- When carrying out oxygen therapy, strictly follow the instructions for use of your O₂ system and your therapy device.

When the pressure measurement ports are not in use, always close them with the cap provided. Otherwise air will escape from the ports and the therapy pressure will not be reached.

3.6 Compatible devices

Therapy devices

JOYCE® masks can be combined with all devices designed for carrying out the therapies described under "Intended Use".

Note that if unsuitable device combinations are used, the actual pressure in the mask may not correspond to the therapy pressure you have been prescribed.

This may also happen even if the correct pressure is indicated on the therapy device. Have your device combination adjusted by a physician or distributor so that the pressure in the mask corresponds to your therapy pressure.

Air humidifiers

JOYCE® masks can be used with either cold-air or warm-air humidifiers.

Never set the humidification level so high that condensation forms in the air hose. Follow the instructions for use of your humidifier and therapy device.

3.7 After use

1. Unhook one of the lower headgear clips from the mask and remove the mask.
2. Loosen the headgear clips and remove the headgear from the mask.
3. Proceed as described in 4 "Cleaning and disinfecting instructions".

3.8 Disassembling/reassembling the mask

Disassembling the mask (Fig. D - A)

1. Detach the mask cushion from the mask frame.
2. Remove the cap from the pressure measurement ports.
3. Rotate the retaining ring counter-clockwise. Now pull the retaining ring and the elbow away from the mask frame.
4. Press the elbow through and out of the retaining ring and pull the rotating sleeve off the elbow.
5. Detach the coarse adjustment component from the mask frame.
6. Detach the fine adjustment component from the coarse adjustment component. Then detach the forehead cushion from the fine adjustment component.

Reassembling the mask (Fig. A - D)

1. Attach the forehead cushion to the fine adjustment component. Then attach the fine adjustment component to the coarse adjustment component.
2. Attach the coarse adjustment component to the mask frame.
3. Push the rotating sleeve onto the elbow until you hear it click into place.
4. Press the elbow into the retaining ring until you hear it click into place.
5. Attach the retaining ring and elbow onto the mask frame. Turn the retaining ring clockwise until you hear it click into place.
6. Push the cap onto the pressure measurement ports.
7. Attach the mask cushion to the mask frame.
8. Fasten the headgear to the mask with the clips.

4 Cleaning and disinfecting instructions

4.1 Schedule

Daily

Clean the mask in warm water to remove all residues. You can disassemble the mask for cleaning. Use a mild cleaning agent (e.g. dishwasher detergent). After cleaning, rinse all parts carefully with clean water.

Caution!

For patients with a weakened immune system or a special medical history, **daily** disinfection or sterilisation may be necessary. In this case, please consult your doctor.

Weekly

Disassemble the mask and clean and sterilise the parts of the mask (see "Cleaning, disinfection, sterilisation" on page 13).

4.2 Cleaning, disinfection, sterilisation

Permissible procedures

Procedure Components	Dishwasher	Washing machine	Disinfection*	Sterilisation (alternative to disinfection)
Mask cushion	•		•	•
Forehead cushion	•		•	
Port cap	•		•	•
Mask frame	•		•	
Coarse adjustment component	•		•	
Fine adjustment component	•		•	
Elbow	•		•	
Rotating sleeve	•		•	
Retaining ring	•		•	
Headgear clips	•	•	•	
Headgear		•		

Caution!

Always clean the components of the mask thoroughly before disinfecting or sterilising them as indicated in the above Table.

Description of procedures

Dishwasher	Rinse cycle at 65 °C. Then rinse all parts carefully with clean water.
Washing machine	Rinse cycle at 30 °C.
Disinfection*	Clean mechanically, e.g. with a brush, in dilute solution. After disinfection rinse all parts thoroughly with distilled water and dry thoroughly.
Sterilisation (alternative to disinfection)	Steam sterilisation in equipment rated according to EN 285. Temperature 134 °C, minimum sterilisation time 3 minutes.

*We recommend the disinfectant GIGASEPT FF. This disinfectant has proven compatibility with the mask. Follow the manufacturer's instructions without fail.

Note

We recommend that you wear gloves (e.g. disposable gloves) while cleaning and disinfecting your mask.

4.3 After cleaning and disinfecting

Dry all parts thoroughly after cleaning and disinfecting. Check the components for residual contamination or damage. If necessary, repeat the cleaning and disinfecting. Replace damaged parts. Reassemble the mask.

4.4 Multiple patient use

If the mask is used by a different patient, it must be cleaned and disinfected first.

Disassemble the mask. Disinfect or sterilize all parts as described in section 4.2 "Cleaning, disinfection, sterilisation". Use a new headgear.

5 Funcional checks

Carry out a visual check after each cleaning and disinfection procedure.

If the mask cushion is damaged and is no longer airtight, replace it with a new one.

Check the functionality and airtightness of the ball-and-socket joint on the elbow. If it no longer moves easily or is no longer airtight, replace the elbow and retaining ring.

Parts that are no longer usable can be disposed of as normal household waste.

6 Troubleshooting

Problem	Cause	Solution
Painful pressure on the face.	Mask is too tight.	Loosen headgear slightly. Adjust position of forehead support.
Draught in the eyes.	Mask is too loose.	Tighten headgear slightly. Adjust position of forehead support.
	Mask does not fit.	Use a different size of mask; contact your distributor.
	High therapy pressure.	Use hard mask cushion (for pressures up to 40 hPa). Available from your distributor.
Therapy pressure is not reached.	Mask not adjusted correctly.	Adjust headgear. Adjust position of forehead support.
	Mask cushion is damaged.	Replace mask cushion.
	Hose system is leaking.	Check push-fit connector and fit of hose sleeves.
	Air is escaping from pressure measurement ports.	Cover pressure measurement ports with cap (included in package).

7 Package contents

Note

All masks are supplied fully assembled (including JOYCE® strap headgear).

JOYCE® models

Article	Article no.
JOYCE®, Size S	WM 26110
JOYCE®, Size M	WM 26120
JOYCE®, Size L	WM 26130
JOYCE® 40 hPa, Size S	WM 26111
JOYCE® 40 hPa, Size M	WM 26121
JOYCE® 40 hPa, Size L	WM 26131

JOYCE® NV models

Article	Article no.
JOYCE® NV, Size S	WM 26160
JOYCE® NV, Size M	WM 26170
JOYCE® NV, Size L	WM 26180
JOYCE® NV 40 hPa, Size S	WM 26161
JOYCE® NV 40 hPa, Size M	WM 26171
JOYCE® NV 40 hPa, Size L	WM 26181

Replacement parts

Article	Article no.
Ball-and-socket joint set: retaining ring, elbow, rotating sleeve	WM 26206
Ball-and-socket joint set NV: retaining ring NV, elbow NV, rotating sleeve NV	WM 26207
Forehead support: Forehead cushion, forehead support coarse adjustment component, forehead support fine adjustment component	WM 26201
Forehead cushion	WM 26200

Article	Article no.
Mask cushion 30 hPa, Size S	WM 26210
Mask cushion 30 hPa, Size M	WM 26220
Mask cushion 30 hPa, Size L	WM 26230
Mask cushion 40 hPa, Size S	WM 26211
Mask cushion 40 hPa, Size M	WM 26221
Mask cushion 40 hPa, Size L	WM 26231
Headgear, blue: JOYCE® strap, incl. headgear clips	WM 26290
Headgear clips (4 pcs.)	WM 26205
Port cap	WM 26208
Instructions for use DE, GB, FR, NL, IT, DK, NO RU (pdf, Internet: www.weinmann.de)	WM 16950 WM 16814
TR (pdf, Internet: www.weinmann.de)	WM 16943

Accessories

Article	Article no.
Headgear, blue: JOYCE® strap EXTRA, 3-point, Size M, incl. headgear clips	WM 26291
JOYCE® strap EXTRA, 4-point, Size M, incl. headgear clips	WM 26292
Head gear, blue: JOYCE® cap, incl. headgear clips	WM 26293
Exhaled air system: Silentflow 2	WM 23600
Noise reduction system	WM 23685

8 Specifications

Technical data

JOYCE®	Value
Device class according to Directive 93/42 ECC:	II a
Dimensions (WxHxD):	80 x 130 x 90 mm
Weight:	c. 110 grams
Therapy pressure range: mask cushion soft: mask cushion hard:	4 to 30 hPa 4 to 40 hPa
Hose connector: cone complying with EN ISO 5356-1 Vented: Non-Vented:	Ø 22 mm (male) Ø 22 mm (female)
Pressure measurement ports:	Ø 4 mm
Width of headgear:	max. 22 mm
Temperature range Operation: Storage:	+ 5 °C to + 40 °C - 20 °C to + 70 °C
Flow resistance Vented: at 50 l/min at 100 l/min Non Vented: at 50 l/min at 100 l/min	0.02 hPa 0.11 hPa 0.09 hPa 0.37 hPa
Sound pressure level at 10 hPa	28 dB(A)
Working life	6 to 12 months*

* depending on frequency of use and cleaning.


Materials

Article	Material
Mask cushion, port cap	Silicone
Forehead cushion, retaining ring	TPE, Polypropylene
Mask frame, elbow	Polycarbonate
Headgear JOYCE® strap	Cotton, polyamide, elasthane
Rotating sleeve, forehead support coarse adjustment component, forehead support fine adjustment component	Polypropylene
Headgear clips	POM

All mask components are latex-free.

Pressure/flow characteristic curve

The pressure/flow characteristic curve shows outlet flow as a function of therapy pressure.

Output pressure of therapy device (hPa)


(Subject to design changes)

9 Warranty

- Weinmann guarantees that the product is free from defects, if used as intended, for a period of 6 months from date of purchase.
- Warranty claims must be accompanied by proof of purchase showing merchant and date of purchase.

- We acknowledge no liability under warranty in the event of:
 - Non-compliance with the instructions for use
 - Incorrect operation
 - Improper use or treatment
 - Repairs to the device attempted by unauthorised persons
 - Force majeure
 - Transport damage due to inadequate packing of returned items
 - Deterioration in normal use and customary wear
 - Failure to use original replacement parts.
- Weinmann does not accept liability for damage consequential upon defects unless due to willful acts or gross negligence or in the event of harm to life or limb due to minor negligence.
- Weinmann reserves the right to choose whether to rectify the defect, to supply an item free of defects or to make an appropriate reduction in the purchase price.
- If a warranty claim is refused, we will not accept the transport cost to or from our address.
- This does not affect your statutory rights.

EEC Letter of Conformity on Medical Devices

We, the manufacturer	Weinmann, Geselle für Medizin GmbH + Co. KG Konsoldweg 40 · D – 22525 Hamburg
declare in sole responsibility the	Product Name: Nasal mask Type/Model: JOYCE® Gr. S / Gr. M / Gr. L JOYCE® NV Gr. S / Gr. M / Gr. L
being in conformity with the respective regulations of the following guideline:	Medical Device Directive 93/42/EEC
Class location:	IIa
Standards used	In particular: DIN EN ISO 17510-2 (May 2003)
Hamburg, 22.07.2005	
	
Dr. C. A. Feldmann Executive Vice President Chief Operating Officer	
© 2004 kb	

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Geräte für Medizin GmbH+Co KG
Geräte für Medizin GmbH+Co. KG