Quick Start Guide

English ..... A

Español ..... B

Français ..... C

# BEFORE YOU START

Please read this guide and the Product Information Manual before use and save both sets of instructions for future reference. The Product Information Manual contains important safety information about the use of this device.


If you feel that your device or any accessories are not operating correctly, please contact your healthcare provider.

**NOTE: This device contains no serviceable parts. Refer all therapy and equipment queries to your healthcare provider.**


## TABLE OF CONTENTS

1. Getting Started .....	A - 2
2. Getting to Know Your F&P ICON™ .....	A - 4
3. F&P ICON™ Menu System .....	A - 6
(1) On/Off and Ramp .....	A - 6
(2) Humidity .....	A - 6
(3) Alarm On/Off .....	A - 6
(4) Alarm Time .....	A - 7
(5) Alarm Volume .....	A - 7
(6) Clock Time .....	A - 7
(7) Brightness .....	A - 7
(8) Sleep Data .....	A - 8
(9) SensAwake™ On/Off .....	A - 8
(10) User Preferences Setting .....	A - 9
(i) ECO Mode .....	A - 9
(ii) Display Pressure .....	A - 9
(iii) Altitude Adjustment .....	A - 9
(iv) Error Code .....	A - 9
4. Communicating Sleep Data .....	A - 10
(1) Telephone Reporting .....	A - 10
(2) Downloading Data to the SmartStick™ .....	A - 10
5. Cleaning and Maintenance .....	A - 10
6. Also Included With Your F&P ICON™ .....	A - 11
7. Replacement Parts .....	A - 11

# 1. GETTING STARTED


Place the device on a level surface, level with or below head height, to prevent water entering the tubing and the device enclosure. Plug the power cord into the power supply.


Remove the Chamber Lid and reach inside the Water Chamber with your fingers to take hold of the circular chamber grip. Pull the Water Chamber up and out of the device.

**CAUTION: Do not fill the Water Chamber while it is in the device; the Water Chamber should always be removed before filling.**


Fill the Water Chamber with distilled water up to the maximum water line, which is marked on the side of the Water Chamber.


**CAUTION: Do not use the device without a Water Chamber.**

To replace the Water Chamber, the arrow on the circular chamber grip should face toward the Elbow.


Fit the Chamber Lid by lining up the large dot on the lid ● with the large dot on the top of the device ●. When the two symbols are aligned, apply downward force and turn clockwise to align the large dot with the small dot on the top of the device.


When the Chamber Lid is locked into place, it should be flush, with no gaps for air to leak. It can now be used as a handle for lifting and carrying.

**NOTE: The Chamber Lid must be correctly locked for optimal pressure delivery.**

The Chamber Lid can be unlocked by twisting in an anti-clockwise direction, then removed by lifting away from the device.


Attach the breathing tube to the Elbow and connect to your mask.

**NOTE: Make sure the electrical connectors on the ThermoSmart™ Breathing Tube come together with the electrical connectors on the Elbow.**

When the selection circle  is at the top of the Menu System, a quick press to the SmartDial™ will start therapy. The Home Screen will light up and the clock-face of the Menu System will rotate to indicate therapy is starting. If Ramp is required, press the SmartDial™ again and hold for 3 seconds to activate the Ramp. The Ramp symbol  will light up on the display.

Turn the SmartDial™ to select and adjust settings. Refer to Section 3 for a description of the Menu System.

**Your device is now ready for use.**

## 2. GETTING TO KNOW YOUR F&P ICON™

### Breathing Tube

All models are compatible with a standard breathing tube, as well as the ThermoSmart™ Breathing Tube.

**WARNING: Do not use the ThermoSmart™ Breathing Tube if there is any visible damage that exposes copper wires.**


### Elbow

The Elbow is where the breathing tube connects to the device. The Elbow can be rotated from side to side, depending on your preferred sleeping position.

### Air Filter

Located on the rear of the device, the Air Filter ensures the air is filtered from particles and dust. To access, first pull out the Filter Grill, then pinch the gauze to remove. The Filter Grill must be pushed back into the cavity to keep the Air Filter in place.


**NOTE: When replacing the Air Filter, the black stripe should face towards the device.**

### USB Port/SmartStick™

The **USB Port** houses the SmartStick™ which records patient compliance, efficacy and flow data. The SmartStick™ must be pushed completely into the USB Port until it clicks into place. To remove, push the end of the SmartStick™ for a spring-release retrieval.

The **SmartStick™** is used to record sleep data from the device memory when it is inserted into the USB Port of the device and to store music for use with AlarmTunes™. Refer to the F&P ICON™ Menu System (Section 3) for more information. The SmartStick™ can be removed and sent to your healthcare provider to view and adjust settings. Changes will be automatically transferred when the SmartStick™ is next inserted into the USB Port of the device.

A line drawing of the device's components. At the top is the Chamber Lid. Below it is the Water Chamber. At the bottom is the SmartDial control panel, which features a rectangular display screen and a circular dial below it. Blue lines with small circles at the end point from the text boxes to the corresponding parts of the device.

### Chamber Lid

The Chamber Lid must be correctly fitted for optimal pressure delivery and can also be used as a handle to carry the device. Twist the Chamber Lid anti-clockwise to unlock and remove.

### Water Chamber

The one-piece Water Chamber holds 420 ml of water and is dishwasher safe. Remove the Water Chamber from the device before filling. Fill to the maximum water line, using distilled water only. Do not use hot water. To place the Water Chamber back into the device, line up the arrow on the circular chamber grip to face toward the Elbow.

### Display Screen/Menu System

The Home Screen is the default display where each point of the Menu System can be viewed and represents the available settings. Once a selection has been made within the Menu System, the Display Screen will dim and return to the Home Screen after 30 seconds. Refer to the F&P ICON™ Menu System (Section 3) for detailed instructions and illustrative examples. If the device is left idle for 10 minutes the Menu System will lock . To reactivate turn the SmartDial™ 180°.


### SmartDial™

Intuitive controls allow one-touch navigation of the Menu System, using a **Press, Turn, Press** principle. For example:

- **Press** the SmartDial™ to turn the device on or off.
- **Turn** the SmartDial™ to find the required setting.
- **Press** the SmartDial™ to enter the setting.
- **Turn** the SmartDial™ to adjust the setting.
- **Press** the SmartDial™ to accept the change.

Refer to the F&P ICON™ Menu System (Section 3) for more detailed instructions and illustrative examples.


### 3. F&P ICON™ MENU SYSTEM

#### On/Off and Ramp

- When the selection circle is at the top of the Menu System a quick press to the SmartDial™ will start therapy.
- To activate Ramp, press the SmartDial™ again and hold for 3 seconds. The Ramp symbol will appear ▽.
- To stop and start pressure delivery, press the SmartDial™ when the selection circle is at the top point.
- The lock symbol 🔒 will appear if the Menu System is left idle for 10 minutes. To reactivate turn the SmartDial™ a half turn in either direction (180°).

**NOTE: In the Auto model if the minimum pressure is set to 4 the Ramp function will be disabled.**


#### Humidity

- To control the level of humidity delivered to the mask, turn the SmartDial™ to the Humidity setting 💧. Press to enter and turn the SmartDial™ again to select the desired humidity level from 0 (off) to 7 (maximum). Press again to accept the change. For example, if the air is feeling too cold increase the Humidity setting in increments of 1 or 2 over 2 to 3 nights until you find the level that is comfortable.

**NOTE: The default setting is Humidity level 4. If using without humidity set the humidity level to zero. The humidity level cannot be adjusted if ECO mode is set to ON.**

- The Boost control offers three additional levels to optimize the Humidity setting. Press and hold the SmartDial™ for 3 seconds while in the Humidity setting 💧 and the Boost symbols 📊 will appear. Turn to select low, medium or high and press again to accept the change.

**CAUTION: The Boost control should only be used in extreme situations when experiencing symptoms from your CPAP therapy. The default Boost setting (medium) will be suitable for most patients.**

#### Alarm On/Off

- Turn the SmartDial™ to the Alarm setting 🔔, press to enter, then turn to select either On ✓ or Off ✗. Press again to accept the change.
- If Alarm On ✓ is chosen, turn the SmartDial™ to select either the default bell 🔔 or the music symbol for AlarmTunes™ 🎵, and press again to accept the change.
- Press the SmartDial™ once to snooze the Alarm for 10 minutes, twice to turn the Alarm off and a third time to discontinue therapy.


**NOTE: The AlarmTunes™ symbol is only visible when the SmartStick™ is inserted. SmartStick™ Studio software is required to utilize the AlarmTunes™ feature.**

- If AlarmTunes™ is selected, turn the SmartDial™ to select the desired music track (from 1 to 5, depending on how many songs are loaded onto the SmartStick™), or if there are two or more songs loaded select Auto for shuffle and press again to accept the change. Refer to the SmartStick™ Studio CD for SmartStick™ Studio software, and instructions for installing and transferring music to your SmartStick™.


#### Alarm Time

- Turn the SmartDial™ to the Alarm Time setting  and press to enter. Turn the SmartDial™ to select the hour, press to accept. Turn to select the minutes and press again to accept the change.


#### Alarm Volume

- To adjust the volume of the alarm, turn the SmartDial™ to the Alarm Volume setting  and press to enter. Turn the SmartDial™ to adjust the volume from quiet to loud and press again to accept the change.
- To adjust the “click” volume of the SmartDial™, press and hold for 3 seconds while in the Alarm Volume setting. The volume symbol  will flash, turn to select the desired volume and press again to accept the change.


#### Clock Time

- Turn the SmartDial™ to the Clock Time setting  and press to enter. Turn the SmartDial™ to select the hour, minutes, 12-hr or 24-hr time. After each selection, press again to accept the change.


#### Brightness

- Turn the SmartDial™ to the Brightness setting  and press to enter. Turn the SmartDial™ to select the Home Screen brightness level and press to accept the change. To adjust the brightness of the display while accessing the Menu System, press and hold the SmartDial™ for 3 seconds and the Active Brightness symbol will appear .

- Turn to select the desired brightness and press again to accept the change.


**NOTE: The Home Screen setting controls the brightness emitted by the LCD when the device is in stand-by mode, and the Active Brightness setting controls the brightness emitted by the LCD whilst navigating the Menu System.**

## Sleep Data

Turn the SmartDial™ to the Sleep Data setting  and press to enter. Turn the SmartDial™ to either the Centers for Medicare and Medicaid (CMS) report  (US specific) or Patient report  and press to enter.

Turn the SmartDial™ to scroll through each of the following statistics:


- Total Compliant Hours *hr*.
- Average Compliant Hours over the last 30 nights *hr*. Press the SmartDial™ to access the average for the last 7 nights and again for the last night.
- Percentage of nights used greater than 4 hours over the last 30 nights  *>4hr*. Press the SmartDial™ to access the percentage for the last 7 nights, and again for the last night.
- The number of nights used since therapy commenced .
- Check Sum .

**NOTE: A  symbol will appear on the Display Screen when the Sleep Data is within the compliance requirements. Once CMS requirements are met, the data within the CMS sub-menu will stop logging.**



- Exit the Sleep Data setting.

**Additional Sleep Data (NOTE: This information may be restricted by your healthcare provider.)**


- AHI for last 30, 7, and last nights.
- Leak for last 30, 7, and last nights.
- 90<sup>th</sup> percentile pressure for last 30, 7, and last nights.
- SensAwake™ for last 30, 7, and last nights.
- Exit the Sleep Data setting.

## SensAwake™ On/Off

- SensAwake™ is activated as a default setting. To deactivate from models with this feature, turn the SmartDial™ to the SensAwake™ setting  and press to enter. Turn the SmartDial™ to select On  or Off  and press again to accept the change.

**NOTE: This setting only applies to the F&P ICON™ Auto.**

## User Preferences


- Turn the SmartDial™ to the User Preferences  setting and press to enter.

### (i) ECO Mode

ECO Mode allows for the device to operate on a reduced power (75 W). Situations for use in this mode might include long-haul air flights and camping.


**NOTE: Humidity output will be significantly reduced in this mode and is not recommended for long-term use.**


- Turn the SmartDial™ to the ECO setting *ECO* and press to enter. Turn to select On  or Off  and press again to accept the change.
- The ECO symbol *ECO* will be displayed continuously on the Home Screen while the Economy Mode is activated.

### (ii) Display Pressure


The Display Pressure setting provides two display options for the Home Screen. Time is the default setting, however this can be changed to show pressure (cmH<sub>2</sub>O or hPa) if required.

- Turn the SmartDial™ to the Display Pressure setting  and press to enter. Turn to select On  or Off  and press again to accept the change.
- If a different unit of measure is required, press the SmartDial™ and hold for 3 seconds. Turn the SmartDial™ to select either cmH<sub>2</sub>O or hPa and press again to accept the change.


### (iii) Altitude Adjustment


- For models without automatic Altitude Adjustment, turn the SmartDial™ to the Altitude setting  and press to enter. Turn the SmartDial™ to select the Altitude level from 0 to 3,000 m and press to accept the change.
- To choose feet, press the SmartDial™ and hold for 3 seconds while at the Altitude Adjustment setting. When the symbol starts to flash, turn to select feet and press to accept the change. Press the SmartDial™ again to enter the menu and turn to adjust to the desired Altitude level. Press again to accept the change.

### (iv) Error Code

If the Error symbol *Error* flashes on the Home Screen, turn the SmartDial™ to the User Preferences  setting and press to enter. Turn the SmartDial™ again to access the Error Code *Error*. Record the number displayed and contact your healthcare provider.

**NOTE: Settings will not change when the device is turned off.**

## 4. COMMUNICATING SLEEP DATA


---

### Telephone Reporting

Many healthcare providers will ask for sleep data over the telephone to track your CPAP therapy progress. Follow the steps outlined under Sleep Data (Section 3, Point 8) which explains how to access this data.

### Downloading Data to the SmartStick™

If your healthcare provider requests a copy of your sleep data from the SmartStick™ please follow the steps below to download this information:

- (1) If the SmartStick™ is already inserted into the USB Port of the device, check that the pressure is off and there is no SmartStick™ symbol  flashing on the Home Screen. The SmartStick™ can now be safely removed by pushing the end for a spring-release retrieval.
- (2) If the SmartStick™ has been removed, ensure the pressure is off before inserting the SmartStick™ back into the USB Port of the device. Wait until the SmartStick™ symbol  has stopped flashing on the Home Screen before removing the SmartStick™ again.
- (3) Data from the SmartStick™ should then be sent in the manner requested by your healthcare provider.

**NOTE: The flashing symbol indicates that data is being transferred from the device memory to the SmartStick™. Removal of the SmartStick™ during this time may cause data corruption.**

## 5. CLEANING AND MAINTENANCE

---

The device should be cleaned as follows:

- Unplug the device from the power supply.
- Wipe the exterior and basin of the device (when the Water Chamber has been removed) with a clean, damp (not wet) cloth and mild dishwashing detergent. Do not use harsh abrasives, solvents or alcohol as these may damage the device.

### Daily

Clean the Water Chamber and the ThermoSmart™ Breathing Tube or standard breathing tube:

- Disconnect the ThermoSmart™ Breathing Tube or standard breathing tube from the mask and the device Elbow.
- Clean either the ThermoSmart™ Breathing Tube or standard breathing tube with warm, soapy water, rinse thoroughly and hang with both ends pointing to the floor to dry.
- Remove the Water Chamber and discard any remaining water.
- The Water Chamber can be washed in a domestic dishwashing machine, or cleaned and washed with mild dishwashing detergent, then rinse and dry.

### Weekly

Thoroughly clean the Water Chamber:

- Soak the inside of the Water Chamber for 10 minutes in a solution of one part white vinegar to two parts water.
- Empty the vinegar solution and rinse the Water Chamber thoroughly with water.

It is recommended to replace your ThermoSmart™ Breathing tube or standard breathing tube, and Water Chamber regularly. Contact your healthcare provider for advice on how frequently your equipment should be replaced. Equipment should be replaced immediately if there is any sign of cracking, discoloration or leaking.

**NOTE: The device is not serviceable and should not be disassembled. Refer all therapy and equipment queries to your healthcare provider.**

## 6. ALSO INCLUDED WITH YOUR F&P ICON™


### F&P ICON™ User Instructions

The Quick Start Guide will explain how to set-up your device for first use, as well as cleaning instructions. The Product Information Manual includes warnings, cautions, product specifications and a trouble shooting guide.

### SmartStick™ Mailer

The SmartStick™ Mailer is used to send your SmartStick™ to your healthcare provider via the post.

### Carry-bag

There are two Carry-bag options for the F&P ICON™. The first is the standard Carry-bag which has been designed for ease of use when packing your CPAP. The second is the Protective Carry-case which provides extra protection when traveling.

## 7. REPLACEMENT PARTS

900ICON200	Water Chamber (10 pack)	900ICON214	Silver Gloss Lid
900ICON202	SmartStick™ (5 pack)	900ICON216	Matt Gray Lid
900ICON203	SmartStick™ Net (5 pack)	900ICON218	Filter Grill
900ICON204	Elbow	900ICON503	Air Filter
900ICON206	Outlet Seal	900ICON308	Protective Carry-case
900ICON208	ThermoSmart™ Breathing Tube*	900ICON315	Carry-bag
900HC221	Breathing Tube*		

\* Applied Parts - to fit 0.86 inches (22 mm) Conical Connector

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Español ..... B

Français ..... C

*This product is covered by one or more of the following patents or patent applications in the corresponding countries listed below.*

**F&P ICON™ Auto** USA: US6598604, US7086399, US6988994, US7186221 Australia: AU 756477, AU2007318331, AU2004266564, AU2006297468 Europe: EP2079504, EP1654024, EP1928529

**F&P ICON™ Premo** USA: US6598604, US7086399 Australia: AU756477, AU2007318331 Europe: EP2079504

**F&P ICON™ Novo** USA: US6598604, US7086399 Australia: AU756477, AU2007318331 Europe: EP2079504

# TABLE OF CONTENTS

---

1. Intended Use .....	A - 1
2. Symbol Definitions .....	A - 1
3. General Warnings and Cautions .....	A - 2
3.1 Contraindications .....	A - 2
3.2 Warnings .....	A - 2
3.3 Cautions .....	A - 4
4. Product Specifications .....	A - 5
5. F&P ICON™ Models and Features Matrix .....	A - 7
6. Device and Consumables Disposal Instructions.....	A - 7
7. Warranty Statement.....	A - 7
8. Troubleshooting .....	A - 8


## 1. INTENDED USE

---

The device is for use on adult patients for the treatment of Obstructive Sleep Apnea (OSA).  
 The device is for use in the home or sleep laboratory.

## 2. SYMBOL DEFINITIONS

---

	Type BF Applied Part		Alternating Current
	ATTENTION Consult accompanying documents		Class II Medical Electrical Equipment
	Do not discard as regular rubbish		CE0123 Conforms with medical device directive 93/42/EEC
<b>IPX2</b>	Drip-proof		Date of manufacture
<b>REF</b>	Catalogue number	<b>SN</b>	Serial number
<b>R</b> Only	Prescription only	<b>ECREP</b>	Authorized representative in the European community
	NZ Radio Interference C-tick mark		Caution
	UL Classified symbol		

## 3. GENERAL WARNINGS AND CAUTIONS

---

This manual refers to the **F&P ICON™** series model as “the device”.

The device treats OSA by delivering a flow of continuous positive airway pressure (CPAP) at a level prescribed by the physician, to splint open the airway and prevent airway collapse.

### 3.1 Contraindications

- Research indicates that the following pre-existing conditions may contraindicate the use of positive pressure for some patients: pneumothorax, bullous lung disease, pneumocephalus, cerebrospinal fluid leak, recent cranial surgery or trauma, abnormalities of the cribriform plate, pathologically low blood pressure or in patients whose upper airways are bypassed.

#### 3.1.1 Precautions:

- The safety and effectiveness of positive pressure has not been established in patients with respiratory failure or chronic obstructive pulmonary disease.
- The safety and effectiveness of the auto-adjusting device has not been established in patients with congestive heart failure, obesity hypoventilation syndrome or central sleep apnea.

#### 3.1.2 Adverse effects:

- Nosebleeds, ear and sinus discomfort may occur from the use of positive pressure therapy.

### 3.2 Warnings

#### 3.2.1 To avoid electric shock from the device:

- Only operate if the device, power cord and plug are dry and in good working order.
- If water damage occurs to your device, disconnect the power cord, discontinue use immediately and seek advice from your healthcare provider.
- Do not store or use the device where it can be pulled into water.

#### 3.2.2 To avoid choking, or inhalation of a foreign object:

- Never place any non-approved objects into any opening of the device, ThermoSmart™ Breathing Tube, Breathing Tube or mask.
- Ensure the recommended filter is fitted to the device before use.
- Ensure the ThermoSmart™ Breathing Tube or Breathing Tube is positioned so it can not become entangled with the body or furniture during sleep.

#### 3.2.3 To avoid burns:

- Do not fill the Water Chamber with hot water.
- After use, wait for the water to cool before touching, carrying or emptying the Water Chamber.
- Do not touch the water in the Water Chamber while the device is operating.
- Position the ThermoSmart™ Breathing Tube so it is uncovered and free from bedding or other materials.

### 3.2.4 To ensure optimal therapy:

- Do not operate the device, Water Chamber, ThermoSmart™ Breathing Tube or Breathing Tube if it is dropped, damaged or not working as intended.
- Pressure adjustments should only be made by a qualified healthcare provider.
- Only use masks, ThermoSmart™ Breathing Tubes, Breathing Tubes and accessories compliant with ISO 17510-2, distributed for use with this device, and recommended by Fisher & Paykel Healthcare or your healthcare provider.
- Do not block the exhaust flow holes of the mask as they are designed to allow a continuous flow of air out of the mask and if blocked CO<sub>2</sub> re-breathing may occur.
- In the event of power failure, machine malfunction or if the device is turned off, remove the mask immediately as the flow through the mask may be insufficient to clear all exhaled gas, which may result in CO<sub>2</sub> re-breathing occurring that can be hazardous.
- On models without automatic altitude adjustment, ensure the altitude level is manually adjusted to ensure optimal pressure delivery.
- Do not use the device without a Water Chamber in place.
- If using the device without water in the Water Chamber it is recommended to set the humidity level to zero.

### 3.2.5 Other:

- Ensure the device is stored and operated on a surface, level or below head height, to prevent water entering the tubing and the device enclosure.
- The device is not intended to be used as a life-supporting device.
- The device complies with the electromagnetic compatibility requirements of IEC60601-1-2. In certain circumstances the device may affect or be affected by nearby portable mobile radio frequency communication equipment, due to the effects of electromagnetic interference. If this should happen, try moving your device or the location of the equipment causing interference, or alternatively consult your healthcare provider.

### 3.2.6 Using supplemental oxygen with your device:

- Supplemental oxygen can be administered at the mask end of the ThermoSmart™ Breathing Tube or Breathing Tube.  
**Note:** At a fixed flow rate of supplemental oxygen the inhaled oxygen concentration will vary, depending on the pressure settings, the patient's breathing pattern, mask selection and leak rate.
- Ensure there is no obstruction downstream of the Oxygen Port as this can affect the delivered oxygen concentration.
- Oxygen concentration should be measured at the point of delivery to the patient.
- Avoid the risk of fire:
  - Only use oxygen when the device is operating. If the device is turned off it can lead to accumulation of oxygen within the device.
  - Ensure adequate ventilation is provided around the device.

- Remove any source of ignition: such as cigarettes, an open flame, or materials which burn or ignite easily at high oxygen concentration.
- Keep oxygen regulators, cylinder valves, tubing, connections and all other oxygen equipment away from oil, grease or greasy substances. Spontaneous and violent ignition may occur if these substances come into contact with oxygen under pressure.

### **3.3 Cautions**

#### **3.3.1 To prevent water damage to the device:**

- Remove the Water Chamber from the device before filling.
- Do not fill the Water Chamber above the maximum level.
- Do not move, carry, transport or store the device with water in the Water Chamber.

#### **3.3.2 General:**

- Only use the device within the Operating Conditions specified in Section 4.
- Position the device so the power cord connection to the power supply is easily accessible.
- Only clean the device in accordance with the cleaning instructions set out in the CLEANING AND MAINTENANCE Section of the Quick Start Guide and only when it is disconnected from the power supply.
- Only use the F&P ICON™ SmartStick™ with the device. Use of any other USB drives may cause data corruption. Do not attempt to change the directories or view the data without software distributed or designed for use with the F&P ICON™.

## 4. PRODUCT SPECIFICATIONS

**DIMENSIONS:** 6.3 H x 6.7 W x 8.7 D in. (160 H x 170 W x 220 D mm)

**WEIGHT:** 4.8 lb (2.2 kg) Packed Weight 8.7 lb (4.0 kg)

**PERFORMANCE: Pressure Range:**  
4 to 20 cmH<sub>2</sub>O (in the unlikely event of fault conditions, pressure may reach up to 29 cmH<sub>2</sub>O)

<b>Maximum Flow Rates</b>					
CPAP Pressure Setting (cmH <sub>2</sub> O)	4	8	12	16	20
Measured flow at patient connection port (Lpm) - Auto/Premo	62.1	101.5	129.4	148.8	143.3
Measured flow at patient connection port (Lpm) - Novo	46.7	52.4	56.5	60.6	64.5

<b>Dynamic Pressure Stability</b>					
CPAP Pressure Setting (cmH <sub>2</sub> O)	4	8	12	16	20
Pressure difference (cmH <sub>2</sub> O) at mask connection port - Auto/Premo	0.52	0.47	0.47	0.52	0.54
Pressure difference (cmH <sub>2</sub> O) at mask connection port - Novo	1.16	1.22	1.28	1.38	1.44

<b>Static Pressure Stability</b>		
	Auto/Premo	Novo
Pressure change (cmH <sub>2</sub> O) at connection port at a pressure setting of 10 cmH <sub>2</sub> O	0.2	0.4

**HUMIDITY:** **With ThermoSmart™ Breathing Tube:**

**Maximum Humidity** = 36 mg/L Body Temperature and Pressure Saturated (BTPS), 82% RH at 10 cmH<sub>2</sub>O, with Humidity level 7 and Boost level high.

**Typical Humidity** = 24 mg/L (BTPS), 90% RH at 10 cmH<sub>2</sub>O, with Humidity level 4 and Boost level medium.

**With Standard Breathing Tube:**

**Maximum Humidity** = 32.24 mg/L, 73.21% RH at 10 cmH<sub>2</sub>O with Humidity level 7 and Boost level high.

**Typical Humidity** = 17.97 mg/L, 85.88% RH at 10 cmH<sub>2</sub>O with Humidity level 4 and Boost level medium.

- GAS TEMPERATURES:** Maximum = 100 °F (38 °C)
- NOISE LEVEL:** Sound Pressure Level < 29 dBA  
Average Sound Power Level < 37 dBA
- WATER CHAMBER VOLUME:** 420 ml up to the maximum fill line
- STANDARDS COMPLIANCE:** IEC60601-1: 1988 + A1 & A2  
AS/NZS 3200.1.0: 1998  
EN60601-1: 1990 + A1-A2  
UL60601-1: 2003

**DATA RECORDING:** The SmartStick™ records up to 365 days of summary efficacy data, 30 days of detailed efficacy data (AHI, Leak), Pressure and Leak data (50 Hz), as well as cumulative data from first use. The device memory can also record up to 365 sessions of summary data or up to 200 sessions (400 hours) of detailed efficacy data.

Electrical Ratings		
Rated Supply Voltage	Rated Current Input	Rated Supply Frequency
110-115 V	1.27 A (1.43 A Max)	50-60 Hz/400 Hz
220-240 V	1.07 A (1.21 A Max)	50-60 Hz

**STORAGE AND TRANSPORT:** The device should always be stored and transported within the following temperatures:

**Transport:** 14 to 140 °F (-10 to 60 °C)

**Storage:** 14 to 140 °F (-10 to 60 °C)

**CAUTION:** Ensure water is emptied from the Water Chamber before traveling. The Protective Carry-case will not close until the Chamber Lid is removed from the device.

- OPERATING CONDITIONS:**
  - Ambient Temperature:** 41 to 95 °F (5 to 35 °C)
  - Humidity:** 5 to 95% Relative Humidity
  - Altitude:** 0 to 9,000 ft (0 to 3,000 m)
  - NOTE:** Above 4,500 ft (1,500 m) the maximum operating pressure will be reduced.

## 5. F&P ICON™ MODELS AND FEATURES MATRIX

Refer to the Quick Start Guide, Section 3, for more information.

Performance Features	Auto		Premo	Novo
	ICONAA <sub>X</sub>	ICONAB <sub>X</sub>	ICONPB <sub>X</sub>	ICONNA <sub>X</sub>
Fully integrated	•	•	•	•
ThermoSmart™ Technology*	•	•	•	•
Auto-Adjusting Pressure	•	•		
Efficacy Reporting	•	•	•	
Compliance Reporting	•	•	•	•
SmartStick™	•	•	•	•
SensAwake™	•			
Proportional Ramp	•	•	•	•
Auto-Altitude Adjusting	Automatic	Automatic	Automatic	Manual
Leak Compensation	•	•	•	
Clock and AlarmTunes™	•	•	•	•
InfoSmart™ Technologies	Advanced	Advanced	Advanced	Basic

**NOTE:** The <sub>X</sub> denotes country variant

\* In some countries the ThermoSmart™ Breathing Tube needs to be purchased as an accessory to activate ThermoSmart™ Technology

## 6. DEVICE AND CONSUMABLES DISPOSAL INSTRUCTIONS

### Device Disposal Instructions

This device contains electronics. Please do not discard as regular rubbish. Dispose of electronics according to local guidelines.

### Consumable Disposal Instructions



Place the mask, breathing tube and Water Chamber in a rubbish bag at the end of use and discard with normal rubbish.

## 7. WARRANTY STATEMENT

Fisher & Paykel Healthcare warrants that the CPAP (excluding consumable items forming part of the CPAP delivery system), when used in accordance with its instructions for use, shall be free from defects in workmanship and materials and will perform in accordance with Fisher & Paykel Healthcare's official published product specifications for a period of 2 years from the date of purchase by the end user. This warranty is subject to the limitations and exceptions set out in detail in <http://www.fphcare.com/osa/cpap-solutions/cpap-warranty.html>

## 8. TROUBLESHOOTING

If you feel that your device is not operating correctly, please refer to the following suggestions. If the problem persists, please consult your healthcare provider. Do not attempt to repair the device yourself.

Problem	Possible Cause	Solution
I cannot turn the pressure on or off	The selection circle may not be at the 12 o'clock position.	Turn the selection circle to the 12 o'clock position and press the SmartDial™ to turn the pressure on or off.
Insufficient air is delivered from the device	The Chamber Lid may not be correctly fitted to the device. Ensure the Water Chamber is in the device.	Refer to Section 1 of the Quick Start Guide.
"Error" is displayed on the Home Screen of the device	An error may have been detected with the device.	Turn the SmartDial™ to the User Preferences  setting and press to enter. Turn the SmartDial™ again to access the Error Code <b>Error</b> . Record the number displayed and contact your healthcare provider for further instructions.
The Protective Carry-case won't close	The Chamber Lid may not be removed from the device. The Elbow may need to be rotated to either side.	The Protective Carry-case has been designed to protect the device during travel. The Chamber Lid must be removed before closing the Protective Carry-case, and water must be removed from the Water Chamber.
AlarmTunes™ does not work	The Alarm may be set to "off". The Alarm may be set to "buzzer". SmartStick™ Studio software may not have been installed. The Alarm volume may not be set loud enough. The SmartStick™ may not be inserted all the way into the USB port.	Refer to Section 3 of the Quick Start Guide for instructions on the Alarm Setting. Once SmartStick™ Studio software has been installed onto your computer and your favorite songs selected, download onto your SmartStick™ and then ensure the SmartStick™ is inserted into the USB Port of the device.
The Heater-Plate and/or Water Chamber base is warm to the touch even though the device isn't being used	The Heater-Plate and Water Chamber are insulated by the device which can cause heat to be retained. This is perfectly normal and should not cause concern.	If concerned about overheating, unplug the device from the power supply when it is not being used.
I think my Water Chamber is leaking because there is a build-up of water on the Heater-Plate	After therapy the device will cool, which may result in condensation forming inside the Water Chamber and pooling on the heater plate.	Before each use remove the Water Chamber and dry the basin of the device with a cloth. If the water build-up becomes excessive, please contact your healthcare provider.
I don't think my humidifier is working	ECO mode may be activated. When ECO Mode is activated the humidity output will be significantly reduced. Check the Humidity setting is at 4 or above. Check your Water Chamber as it may be empty.	To disable the ECO Mode, turn the dial to the User Preferences  Setting and press to enter. Turn to select <b>ECO</b> and press to enter. Turn to select <b>OFF</b>  and press to accept the change. Refer to Section 3 of the Quick Start Guide for more information on the humidity controls.
The SmartStick™ symbol  is flashing on the Home Screen	The SmartStick™ may not be inserted correctly.	Take the SmartStick™ completely out of the USB port, then insert again until it clicks.
I can't select anything on the Menu System	The screen lock  may have activated after 10 minutes of therapy or no-use.	Turn the SmartDial™ 180° to unlock the Menu System.
The device seems noisy; there is a hissing noise coming from the device	The Chamber Lid may not be properly fitted, causing air to leak.	Remove the Chamber Lid and replace again so that it is completely flush with the device, by following the instructions set out in Section 1.

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