DeVilbiss
HEALTHCARE

DV54 Series

CE 0044

EN

DeVilbiss® SleepCube® Positive Airway Pressure Device

CAUTION—Federal (U.S.A.) law restricts this device to sale by, or on the order of a physician.

Assembled in USA of US & Imported Parts. Instruction Guide.

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TABLE OF CONTENTS

Symbol Definitions	3
Important Safeguards	3
Introduction	
Intended Use/Indications for Use	5
Contraindications.....	5
What Is Obstructive Sleep Apnea (OSA).....	5
How CPAP Therapy Works.....	6
Acclimating To Your Therapy.....	7
Benefits Of Your Therapy	7
Helpful Resources	8
Key Features	
Unpacking The Contents	9
SleepCube Device.....	9
Keypad	9
System Assembly without Humidification.....	10
Keypad and Display.....	10
Operation	
Start Up.....	11
Using The Comfort Delay Feature	11
Using The SmartFlex Feature	12
Shut Down	13
Advanced Menu Controls	13
Advanced Menu List.....	14
SmartFlex Menu List	16
SleepCube Enable Menu List	17
Blower Off Menu	18
SmartCode and Adherence Score Information.....	18
AutoAdjust Patient Menu	19
CPAP Patient Menu.....	19
Quick View Menu Layout	20
Patient Messages	21
Reminders	21
Notifications	21
Travel Information	22
SleepCube DC Operation.....	22
Battery Power	22
Supplemental Oxygen	23
Accessories/Replacement Items	24
Maintenance	25
Cleaning	26

Troubleshooting.....	27
Specifications	29
DeVilbiss Guidance and Manufacturer's Declaration.....	31
Declaration of Conformity.....	35

SYMBOL DEFINITIONS

	Class II electrical protection– double insulated		Type BF equipment– applied part
	Data Port Input/Output		Standby – turns the blower on or off
	Next Item on LCD display		Increase value shown on LCD display
	Previous Item on LCD display		Decrease value shown on LCD display
	Attention - Consult Accompanying Documents		Heat
	DC Input: 12 VDC, center pin is positive		100 – 240 Volts AC input 50/60 Hz
	Locked		Unlocked
	IPX1 Drip-Proof, Vertical		
	Delay – starts the delay cycle, if prescribed		
	This device contains electrical and/or electronic equipment that must be recycled per EC Directive 2002/96/EC - Waste Electrical and Electronic Equipment (WEEE)		

IMPORTANT SAFEGUARDS

When using electrical products, basic safety precautions should always be followed. Read all instructions before using this device. Important information is highlighted by the following terms.

DANGER Urgent safety information for hazards that will cause serious injury or death.

WARNING Important safety information for hazards that might cause serious injury.

CAUTION Information for preventing damage to the product.

NOTE Information to which you should pay special attention.

**PLEASE READ ALL INSTRUCTIONS BEFORE USING
THIS DEVICE.**

SAVE THESE INSTRUCTIONS

DANGER

Electric Shock Hazard – Do not use while bathing.

Electric Shock Hazard – Do not immerse this device into water or any other liquid.

Electric Shock Hazard – Do not attempt to open or remove the enclosure. There are no user-serviceable internal components. If service is required, return the product to your home care provider. Opening or tampering with the product will void the warranty.

WARNING

The DeVilbiss SleepCube® should be used only with masks recommended by DeVilbiss, your physician or respiratory therapist.

To avoid rebreathing of exhaled air, do not use a CPAP mask unless the device is turned on and providing a supply of air. Venting in the mask should never be blocked. When the device is turned on and providing a fresh supply of air, exhaled air is flushed out of the mask vent. However, when the device is not operating, exhaled air may be rebreathed. Rebreathing of exhaled air for longer than several minutes can in some circumstances lead to suffocation. This warning applies to most CPAP devices.

The DeVilbiss SleepCube is not a life support device and may stop operating with certain device faults or with a power failure. It is intended to be used on spontaneously breathing individuals weighing 66 lbs/30 Kg or greater.

To avoid electric shock, always unplug power cord from wall outlet power source when performing cleaning.

Use only accessories recommended by DeVilbiss.

CAUTION–*The circular data port connector located on the back of the SleepCube is used to attach accessories to the device. The connector must only be used with accessories approved for use by DeVilbiss. Do not attempt to attach any other device to this connector as it may damage the CPAP or the accessory device.*

INTRODUCTION

CAUTION—*Never rinse or place the device in water. Never allow liquids to get into or around any of the ports, switches or air filter; doing so will result in device damage. If this occurs, discontinue use and remove the power cord from the power source. Allow the device to completely dry before use.*

CAUTION—*Do not place the SleepCube where it can be bumped onto the floor or where the power cord may create a trip hazard.*

CAUTION—*Only the DeVilbiss DV5 series Heated Humidifier system is recommended for use with the SleepCube. Other humidifier systems may prevent the device from detecting snoring and may cause inappropriate pressure levels in the mask.*

CAUTION—*Oxygen is a prescription gas and should only be administered under the supervision of a physician.*

INTRODUCTION

Intended Use/Indications for Use

The DeVilbiss® SleepCube® Model DV5 Series is intended for use in treating obstructive sleep apnea in patients 66 lbs/30 kg and above.

Contraindication

Do not use SleepCube if your upper airway has been bypassed.

NOTE—*Symptoms of dryness of the throat, nasal passages, and the mouth are common with positive airway pressure therapy. This is especially true in dry climates and during the cold season when humidity in the air is typically lower than at other times. The SleepCube features an optional humidifier system to help minimize these effects.*


What Is Obstructive Sleep Apnea (OSA)?

Obstructive Sleep Apnea (OSA) is a breathing disorder that affects more than 22 million Americans – approximately 4% of middle-aged men and 2% of middle-aged women.

Physiological factors such as size of tonsils, shape of the nose, narrowing of the upper airway, or any combination of these can contribute to the potential for having OSA. During sleep, your airway can relax so much that it begins to collapse and becomes obstructed. Relaxation allows the tissue at the back of your throat to vibrate as you breathe causing snoring. But as your airway collapses more fully, it blocks your airflow and stops your breathing altogether. (In fact, the word “apnea” is derived from a Greek term meaning “without breath.”)

This stoppage occurs repeatedly throughout the night, each episode lasting from just a few seconds to more than a minute. As your blood oxygen level drops to dangerous levels, your

central nervous system triggers a sudden gasp for breath, partially awakening you and preventing you from experiencing the sustained, deep sleep your body needs. As a result, it is common to feel fatigued, stressed and irritable the following day – a less-than-optimal and potentially dangerous condition in which to work, drive, and make decisions.

There are important medical concerns as well. OSA starves your heart, brain, and organs of life-sustaining oxygen. Over time, this deprivation can seriously impact your health potentially resulting in memory loss, impotence, hypertension, coronary disease, strokes, and heart attacks. Recent studies indicate that OSA has the same risk factor for heart disease as smoking, high cholesterol and alcohol, and estimate that as many as 38,000 people a year die from the effects of OSA.

Clearly, OSA is a serious medical condition that requires attention and treatment.

How CPAP Therapy Works

There are several treatments for OSA including surgery and CPAP therapy. However, CPAP therapy has proven very effective for most patients. CPAP stands for “continuous positive airway pressure” which is exactly what the DeVilbiss SleepCube device delivers to your upper airway, acting as an “air splint” to keep the passage open while you sleep.

This reduces or eliminates the obstruction, allowing you to enjoy the deep, refreshing, uninterrupted sleep you need. More importantly, your body will get the rest and oxygen it needs with more continuous, uninterrupted sleep – helping you avoid some of the serious health risks associated with OSA. And the SleepCube is designed to provide this therapy to you quietly and comfortably so that you can sleep through the night without distraction.

The DeVilbiss SleepCube AutoAdjust® system automatically adjusts air pressure by measuring your snoring and breathing patterns on a breath-by-breath basis, and

compensating to meet your pressure needs. Your breathing patterns are continuously analyzed, stored and used to make any needed pressure adjustment.

The rate of pressure adjustment can be as high as 10 cmH₂O in a 10 minute period with changes made gradually enough to prevent disturbing your sleep yet quickly enough to respond to your pressure needs. In the absence of respiratory events, the pressure is reduced to maximize your comfort. The SleepCube AutoAdjust is also capable of operating in CPAP mode – without the automatic functionality described above – to give your physician maximum versatility in prescribing treatment for your sleep condition.

Acclimating To Your Therapy

Naturally anything new – including CPAP therapy – takes some getting used to. Having your mouth or nose (or both) covered by a mask may make you feel anxious or claustrophobic. Feeling the rush of air from your SleepCube may cause discomfort or concern. But don't worry – every successful, long-term CPAP user began where you are today with the same reactions and concerns, and overcame them to enjoy the full benefits of their therapy.

Here are some steps to help you gradually ease into your therapy and overcome your anxieties:

1. Try attaching the mask to the CPAP unit, and switching the unit “On.” Practice breathing through the mask for one hour while watching television, reading, or performing some other sitting activity.
2. If you nap, use the CPAP during this time.
3. Use CPAP at least during the first initial 3-4 hours of nighttime sleep.
4. Use CPAP through an entire night of sleep.

Advance through the above steps one step at a time and progress to the next step as soon as you can carry out the current step without anxiety.

Benefits Of Your Therapy

As a DeVilbiss SleepCube user, you are on your way to living a happier, healthier life – a life free from the sleep disruptions and health risks associated with OSA. After you become acclimated to CPAP therapy, you may be surprised at how much better you feel, both physically and emotionally. CPAP therapy may give you the opportunity to finally experience the deep, restful sleep you've been missing out on. Many patients find they have more energy, better judgment, and improved concentration during the day. With effective therapy you'll be less stressed, less irritable and more like your old self again, with a healthy, positive outlook on life. In fact, you may well wish you'd gotten your SleepCube sooner!

Helpful Resources

As you begin your SleepCube therapy, remember: you're not in this alone! You're part of a large and growing network of CPAP users, and many resources are available to answer your questions, address your concerns, and help make your therapy as effective and effortless as possible.

American Academy of Sleep Medicine

One Westbrook Corporate Ctr, Ste 920
Westchester IL 60154
Phone: (708) 492-0390
Fax: (708) 492-0943
www.aasmnet.org

American College of Chest Physicians

3300 Dundee Rd
Northbrook, IL 60062-2348
Phone: (847) 498-1400 or (800) 343-2227
Fax: (847) 498-5460
www.chestnet.org

The Sleep Apnoea Trust

12a Bakers Piece,
Kingston Blount
Oxon.
OX39 4SW
Tel: +44 845 60 60 685
www.sleepapnoea-trust.org

National Heart, Lung & Blood Institute Information Center

P.O. Box 30105
Bethesda, MD 20824-0105
Phone: (301) 592-8573
Fax: (240) 629-3246
www.nhlbi.nih.gov

American Academy of Neurology

1080 Montreal Avenue
St. Paul, Minnesota 55116
Phone: (651) 695-2717 or (800) 879-1960
Fax: (651) 695-2791
www.aan.com

American Sleep Apnea Association

1424 K Street NW, Suite 302
Washington, DC 20005
Phone: (202) 293-3650
Fax: (202) 293-3656
www.sleepapnea.org

National Center on Sleep Disorder Research National Heart Lung & Blood Institute (NIH)

6701 Rockledge Drive
Bethesda, MD 20892
Phone: (301) 435-0199
Fax: (301) 480-3451
www.nhlbi.nih.gov/sleep

National Sleep Foundation

1522 K Street, NW, Suite 500
Washington, DC 20005
Phone: (202) 347-3471
Fax: (202) 347-3472
www.sleepfoundation.org

KEY FEATURES

KEY FEATURES

Unpacking The Contents

The following items are included in your DeVilbiss SleepCube system:

- | | | |
|----------------------|----------------------|--------|
| 1. SleepCube device | 4. Carrying case | 7. DVD |
| 2. Spare filter pack | 5. Air supply tubing | |
| 3. AC line cord | 6. User manual | |



If any of the above items are missing, please contact your home care provider.

NOTE—If you received a humidifier with your SleepCube system, please refer to the humidifier's instruction guide for additional information.

DeVilbiss SleepCube Device

(Figures A and B)


1. Keypad/LCD Display (see Figure C)
2. Air Supply Port on Back
3. Air Supply Port on Bottom (for optional humidifier)
4. Air Supply Port Plug
5. Heater Power Connector (for optional humidifier)
6. AC Power Connector
7. DC Power Connector
8. Data Port
9. Power Cord (appropriate to your wall outlet) (not shown)
10. Air Inlet Filter Opening
11. Connector Cover (not shown)

Keypad (Figure C)

- | | |
|------------------|---|
| 1. On/Off | 5. Decrease Value |
| 2. Previous Item | 6. Increase Value |
| 3. Next Item | 7. Heater Power LED (for optional humidifier) |
| 4. Delay | |

Ask your home care provider for information about other DeVilbiss masks, equipment, and accessories.

SYSTEM ASSEMBLY WITHOUT HUMIDIFICATION

NOTE—If your system includes a humidifier, please refer to the instructions provided with the humidifier for assembly and operation.

1. Locate the air supply port and heater accessory connector on the bottom of the SleepCube device. Ensure that each opening has the appropriate soft cover securely attached.
2. Place the SleepCube on a stable surface such as a nightstand or table. If you wish, you may also place the device on the floor at the bedside, being careful to place the unit where it won't be kicked or stepped on. Ensure that the air inlet in the back of the device is not blocked by anything such as curtains or bedding.


NOTE—Never place the SleepCube system on a soft surface such as a bed or couch during operation.

CAUTION—The height of the SleepCube system must always be lower than the interface mask when using a humidifier to prevent water from getting into the CPAP mask.

3. Insert the power cord into the back of the device. Make sure the power cord is fully inserted into the connector. Plug the other end of the power cord into the AC wall outlet.
4. Attach the air supply tubing to the air supply interface in the back of the device.
5. Prepare the mask (supplied by your home care provider) by referring to the mask's instruction guide.
6. Connect the air supply tubing to the mask.

NOTE—The air supply tubing supplied with the device is specially designed to prevent kinking and allow you to receive optimal airflow. Use only the 22 mm diameter, smooth-bore, kink-resistant tubing supplied with this device or an air supply tubing recommended by your home care provider.

KEYPAD AND DISPLAY


Basic DeVilbiss SleepCube Controls		
Name	Symbol	Function
ON/OFF		The ON/OFF button starts and stops the flow of air from the device.
HEATER		If your SleepCube has a heater/ humidifier attached, these buttons adjust the heater setting. Increase the setting for more humidity and decrease the setting for less humidity.
DELAY		The Delay button activates the comfort delay feature and only functions if this feature has been enabled. To understand how to set the parameters of this feature, please refer to the Using the Comfort Delay Feature found on page 11.

NOTE—The additional control buttons that you see on the SleepCube control panel provide advanced controls for specific settings. Please refer to the ‘Advanced Menu Controls’ section on page 13 for more details.

OPERATION

Startup


To begin operation, verify that your SleepCube is connected to a ‘live’ power source. When plugged in, the LCD display should appear as noted.

Put on your CPAP mask and press the ON/OFF button on the keypad to turn on airflow OR put on your CPAP mask and take a few breaths to turn on airflow automatically. The pressure in your mask will rise to the set pressure in about 10 seconds.

NOTE—If the device does not automatically start after taking 2-3 breaths, this feature may have been disabled. If Auto ON/OFF is disabled, you must manually press the ON/OFF button to start and stop the device.

The keypad display shows the actual mask pressure while the device is in use. Because this is the actual pressure, the device pressure reading may vary slightly as you breathe. Your display should appear as noted.

If your DeVilbiss SleepCube system includes a humidifier, its heater controls and status are displayed as noted in the image.


NOTE—If the display on your device is not similar to those listed above, refer to the chapter section titled “Patient Messages” on page 21.

Using The Comfort Delay Feature

Your SleepCube has a comfort delay feature to help make it easier for you to fall asleep. The airflow will start at a low pressure and then gradually ramp-up to the prescribed setting over the specified delay time.

NOTE—When the delay comfort feature is activated, a message will flash on the display periodically, indicating the amount of time remaining in the delay.

You can adjust the amount of time spent in comfort delay by using the advanced device controls. The Delay Time can be adjusted from 0 to 45 minutes in 5-minute increments. Refer to page 14 for details.

NOTE—If you are unable to change the delay time, your provider has locked this setting. Refer to the keypad menu guide on page 14 for details.

Comfort Delay Feature In Standard CPAP Mode

When you use the comfort delay feature while your device is operating in CPAP mode, the airflow you receive at the beginning of your therapy is temporarily lowered to a pressure selected by your physician when the device was setup.

Press the delay button to activate the comfort delay feature. The airflow pressure will drop to the prescribed comfort delay pressure.

During the delay period, the pressure will gently rise to the pressure prescribed for your therapy. You may repeat the comfort delay by pressing the delay button again if you have not fallen asleep.

Comfort Delay Feature in AutoAdjust Mode

When you use the comfort delay feature while your SleepCube is operating in AutoAdjust mode, the airflow you receive at the beginning of your therapy is temporarily lowered to a pressure selected by your physician when the device was setup.

Press the delay button to activate the comfort delay feature. The airflow pressure will drop to the prescribed comfort delay pressure.

During the delay period, the pressure will gently rise to the pressure setting prescribed for your therapy. You may repeat the comfort delay by pressing the delay button again if you have not fallen asleep. The AutoAdjust will not adjust the pressure setting during the comfort delay period, allowing you to fall asleep without responding to early respiratory events.

Using The SmartFlex Feature

Your SleepCube has a SmartFlex pressure relief feature to help make it easier for exhale against your prescribed pressure by lowering the pressure slightly when you exhale. The pressure returns to your prescribed pressure during inhale to maintain your therapy.

P:10.0	FLX	Auto
	Heat:Off	

SmartFlex pressure relief includes three settings for exhale relief; 1, 2, and 3. (1 is the least amount of relief, 3 is the most amount of relief). Setting SmartFlex pressure relief to zero (0) turns the SmartFlex feature OFF.

The pressure waveform can also be adjusted independently for inhale and exhale to smooth the transition when SmartFlex pressure relief is active. The pressure waveform settings are SmartFlex-i (IRnd) and SmartFlex-e (ERnd) with settings from 0 to 5 in steps of 1 (0 is the least amount of smoothing, 5 is the most amount of smoothing).

SmartFlex pressure relief can be operated Full Time, during the Delay Only or turned off.

The display will show “FLX” in the top row of the default screen (between CPAP pressure and Mode) while SmartFlex pressure relief is operating.

When SmartFlex pressure relief is turned on Full Time, the pressure relief operates at all times when you are breathing on the device. When SmartFlex pressure relief is set to Delay Only, the pressure relief operates only during the Comfort Delay period and turns off after the delay period is finished. When SmartFlex pressure relief is turned OFF, the pressure relief does not operate.

The device contains a SmartFlex disengagement feature. If a patient is experiencing a respiratory event while the CPAP is in the relief pressure mode, after a limited time the CPAP will shift to the prescription pressure and remain at this pressure until normal breathing resumes.

Shut Down

Press the ON/OFF button to turn the device off. You may also simply remove your mask. After a few seconds the device will display the message “Mask Off. Please check mask fit.” If you do not put the mask on again the device will automatically turn off in about 20 seconds. Pressing any button will remove the message from the display.

NOTE—If the DeVilbiss SleepCube does not turn off automatically, the “Auto-OFF” feature may have been disabled. In this case you must manually turn the device off.

ADVANCED MENU CONTROLS

With the device on, press the “Next” or “Previous” buttons on the keypad to scroll through the items on the display. Not all of the features in the table below are available in all models; some are optional and some depend on what your physician and/or home care provider has selected. The table shows the items you may see, and describes the feature associated with each item. Also, some features are only available in AutoAdjust Mode or with the DeVilbiss SmartLink™ Therapy Management Module – these are noted in the ‘Mode’ column below. If you do not see some of these options, your model or configuration does not support them.

Name	Symbol	Function
PREVIOUS ITEM		Pressing this button displays the previous menu option.
NEXT ITEM		Pressing this button displays the next menu option.
SETTING CHANGE		These buttons allow you to adjust the setting that is shown on the display.

Please note that the values displayed are examples only.

Advanced Menu List

Display Example	Mode	Description
	All	This is the default screen when the SleepCube is turned on. The actual mask pressure and heater setting are shown. The displayed pressure may vary slightly while breathing. If the heater system is not connected, the heat setting will be blank.
	AutoAdjust Mode	This is a prescription setting which limits the maximum pressure at which the SleepCube will operate if in AutoAdjust mode. Your home care provider programs this setting for you during device setup.
	AutoAdjust Mode	This is a prescription setting which sets the minimum pressure at which the SleepCube will operate in AutoAdjust mode. Your home care provider programs this setting for you during device setup.
	CPAP Mode	This is a prescription setting which controls the operating pressure.
	All	Selectable delay time setting. Your provider has set the delay time.
	CPAP Mode	This is a prescription setting which controls the operating pressure during the delay time.
	All	The standard tubing length supplied with the device is 6 feet (1.8 meters). However, you can adjust your SleepCube to accommodate an air supply tubing length of 10 feet (3.0 meters). Replacement tubing can be obtained through your home care provider (part numbers listed at the end of this manual).
	All	Pressing the down key enters the SmartFlex menu. Refer to SmartFlex Menu section for details.
	All	Pressing the down key enters the Enable menu. Refer to next section for Enable Menu details.
	All	This setting changes the contrast for the LCD display. This only needs to be adjusted if you find it difficult to see the display from your viewing angle.

Display Example	Mode	Description
Compliance Meter 123456.7 hours	All	The total number of hours (in 1/10 hour increments) that the SleepCube has been providing therapy.
Firmware Version V0.20 30/06/2010	All	Shows version of installed firmware.
Serial Number HDI23456	All	Shows the serial number of the device.
Adherence Score 83% (25/30)	All	Shows the maximum adherence achieved for any period of thirty consecutive days of use within a ninety-day time frame. The score represents the percentage of days, within thirty consecutive days, where usage met or exceeded the 'Usage Threshold.'
SmartCode 1-Day F6CD-37F-357R SmartCode 7-Day 3CCC-377D-7J7X SmartCode 30-Day DWY4-3DB7-DJ7L SmartCode 90-Day 752TT-47DB-7DT7C SmartCode Usage ZDN-79FK-K7H	All	SmartCode™ data is used for therapy management reporting; therapy information is encoded in a series of letters and numbers, and includes a special character for error checking to ensure accuracy. Use the UP arrow or DOWN arrow keys to change the SmartCode reporting period. There are five levels of SmartCode data – 1-day, 7-day, 30-day, 90-day and Usage. Use UP or DOWN button to change between levels.

Note- If the unlocked symbol is shown, you are able to adjust this setting. If the locked symbol is shown, this setting can only be adjusted by your provider.

SmartFlex Menu List

The SmartFlex menu is available in all modes and is a sub-menu that allows you to control various functions of your SleepCube. It is accessed by pressing the down button while in this Advanced Menu.

Display Example	Description
	This setting controls the amount of pressure relief. OFF = SmartFlex feature is OFF 1 = least amount of pressure relief 3 = most amount of pressure relief
	This setting controls the SmartFlex operating mode. Full Time means SmartFlex feature is always active Delay Only means SmartFlex feature is active only during the pressure delay cycle.
	This setting controls the Inhale pressure waveform.
	This setting controls the Exhale pressure waveform.
	Pressing the up key exits the SmartFlex menu

Note- *If the unlocked symbol is shown, you are able to adjust this setting. If the locked symbol is shown, this setting can only be adjusted by your provider.*

SleepCube Enable Menu List

The Enable menu is available in all modes and is a sub-menu that allows you to control various functions of your SleepCube. It is accessed by pressing the down button while in this Advanced Menu.

Display Example	Description
	This option controls the Auto-OFF feature, which automatically stops the flow of air after taking off the mask. The device will require at least 30 seconds to turn OFF depending on the kind of mask used. When Auto-OFF is enabled, Auto-ON is also enabled.
	This option controls the Auto-ON feature, which automatically starts the flow of air after breathing once or twice into the mask. This is always enabled when Auto-OFF is enabled.
	This option controls the mask fit check feature. The air flow to your mask is constantly measured and if the amount of airflow exceeds 95 liters/minute for more than 10% of the time used, a notification will be created. This notification will be displayed the next time you use the SleepCube.
	This is a setting which controls the backlight. Low backlight Enable will keep the backlight dim during operation. Low backlight Disabled will turn the backlight OFF during operation.
	Pressing the up key exits the Enable menu.

Note- If the unlocked symbol is shown, you are able to adjust this setting. If the locked symbol is shown, this setting can only be adjusted by your provider.

Blower Off Menu

Compliance information can be viewed without turning the CPAP on. With the blower OFF (display shows OFF), press the left or right arrow keys on the keypad to scroll through the information.

Display Example	Description
OFF	Blower OFF, the CPAP unit is plugged in to AC or DC power with the blower OFF.
Quick View Menu Press ▼	Pressing the down key enters the Quick View menu. Refer to page 20 for Quick View Menu details. Note: <i>The "Quick View" menu is only shown if enabled from the clinical menu.</i>
Serial Number HD000005	Shows the serial number of the device.
Adherence Score 83% (25/30)	Shows the maximum adherence achieved for any period of thirty consecutive days of use within a ninety-day time frame. The score represents the percentage of days, within thirty consecutive days, where usage met or exceeded the 'Usage Threshold.'
SmartCode 1-Day F6CD-37F-357R	SmartCode™ data is used for therapy management reporting; therapy information is encoded in a series of letters and numbers, and includes a special character for error checking to ensure accuracy. Use the UP arrow or DOWN arrow keys to change the SmartCode data reporting period. There are five levels of SmartCode – 1-day, 7-day, 30-day, 90-day and Usage. Use UP or DOWN button to change between levels.
SmartCode 7-Day 3CCC-377D-7J7X	
SmartCode 30-Day DWY4-3DB7-DJ7L	
SmartCode 90-Day 752TT-47DB-7DT7C	
SmartCode Usage ZDN-79FK-K7H	


SmartCode Data and Adherence Score Information

Your Healthcare provider may contact you to retrieve SmartCode data and Adherence Score information. To display the SmartCode data, press the LEFT arrow key once. Use the UP arrow or DOWN arrow keys to change the SmartCode reporting period. Press the LEFT arrow key again to display the Adherence Score.


DV54 AutoAdjust Patient Menu Layout

The following chart shows by Mode each Display you will scroll through using the Next or Previous buttons on the keypad. Pressing the Down key at the Enable Menu Display will allow you to enter this submenu - then use Next or Previous buttons to scroll through it; press the Up button to return to main menu. Once the final Display is reached, pressing the Next button will return you to the first Display.

MENU - AutoAdjust Mode

MENU - CPAP Mode


DV54 Quick View Menu Layout

The Quick View Menu shows a summary of information which is reported in the SmartCode data. Information in the Quick View includes: Number of days used more than 4 hours (or can be set to 5 hours by clinician), AHI (Apnea/Hypopnea Index), High Leak percent (percent of time at leak >95 L/m), and Daily Usage (average hours per day used). Each parameter can be viewed in any of four time periods, 1 day, 7 days, 30 days and 90 days.

The following chart shows the Quick View information that is available on your CPAP display if viewing this information has been enabled by your provider. The Quick View menu shows a summary of your SmartCode data for 1, 7, 30, and 90 day periods. Each of the time periods contains info for “Number of Days > x hrs”, Apnea Hypopnea Index (AHI), High Leak % (>95 L/m) and Average Daily Hours of while breathing use.

Diagram of CPAP displays

PATIENT MESSAGES

Your SleepCube will alert you of issues that may require some action on your part in order to make your therapy more effective. These messages are displayed when the device is turned on, either by pressing the ON/OFF button or by breathing into the mask. You can press any key or button to clear the message from the display, but the SleepCube will operate normally even if the message is not cleared.

Two kinds of messages are displayed: Reminders and Notifications.

Reminders

Reminders alert you when parts of your system need to be replaced. Some components of your system wear over time and, if not replaced, may compromise your therapy. If you see a message on the display that is not listed below, you may have the optional DeVilbiss SmartLink™ Module attached to your device, which provides additional messages. Refer to the documentation that was provided with the SmartLink Module.

Clean Filter – The filter should be checked every 10 days for signs of dirt or wear so that it can be cleaned as needed. This message helps to remind you to check it regularly.

Notifications

Notifications identify conditions in your SleepCube that may require action by you or your home care provider in order to maintain a high level of therapy. If you see a message on the display that is not listed below, you may have the optional DeVilbiss SmartLink Module attached to your device, which provides additional messages. Refer to the documentation that was provided with the SmartLink Module.

NOTE—*The values displayed are examples only.*

Delay Running xx Minutes Left – While a comfort delay session is active, this message is flashed on the display every 5 seconds to let you know how much time is left in the delay.

Mask Leak – This message means the SleepCube has detected a large amount of airflow during therapy for at least 10% of the time during your previous use session. This notification is displayed when the device is turned on. If this message is displayed, put on the mask and adjust the headgear to ensure the mask is properly fitted to your face. Follow the mask manufacturer's directions for adjusting the mask and headgear for proper fitting. Available on AutoAdjust models.

Mask Off – This message appears whenever a high airflow condition occurs. High airflow is normally due to a poor mask fit

or a removed mask. Check for leak around the mask seal and make adjustments as necessary according to the mask manufacturer's instructions. This message will be displayed until the high airflow problem is corrected. If the high airflow condition persists for about 20 seconds, the device will automatically turn off (if AutoOFF is enabled).

Device Fault – If a device fault message is displayed, refer to Troubleshooting for instructions.

Device Fault	E07
Call Provider	

SLEEP CUBE TRAVEL INFORMATION

International Power Changes

Your DeVilbiss SleepCube is equipped with a universal power supply. This means that it is automatically capable of accepting line voltages of 100 to 240V~ 50/60Hz. You do not need to make any adjustments to the device if you are traveling to another country where the line voltage is different than your domestic line voltage. Simply contact your home care provider for the correct style power cord for the region in which you will be traveling (refer to Accessories/Replacement Items for the correct part number).

Altitude

Your SleepCube automatically compensates for changes in altitude between sea level and 9000 ft (2750 m). You do not need to make adjustments to the device for changes in altitude.

SLEEP CUBE DC OPERATION

Battery Power

Your SleepCube can be operated from a 12V DC power source in a motor home, boat or other recreational vehicle. The power source can be used with a 12V DC accessory cable or a 12 V DC stand-alone battery. Both parts are available from your home care provider (refer to Accessories/Replacement Items for the correct part numbers). Both accessory cables are required for stand-alone battery operation. The stand-alone battery connection uses clamps that attach to the terminal connections of the battery.

NOTE—*The heater for the SleepCube's humidification system will not operate directly from a 12V DC power source. If you need humidification, non-heated humidification is still available. Alternately, instead of connecting 12V DC power directly to the SleepCube, you can pass the DC power through an inverter that will deliver AC power to the system. This will allow the heater to function. The power rating of the inverter must be at least 200 watts @ 110V~ / 400 watts @ 220V~.*

NOTE—*The SleepCube will automatically switch power sources between AC and 12V DC if*

both sources are available. If your unit is connected to both AC and 12V DC power sources, it will operate on AC power as long as AC power is present, and automatically switch to 12V DC power if AC power is lost. When AC power is restored, it will return to the AC source. You do not need to reset or adjust any settings when the power source switches, but remember that your heater will not operate on 12V DC.

SUPPLEMENTAL OXYGEN

Supplemental oxygen is sometimes prescribed in addition to CPAP for treatment of Obstructive Sleep Apnea. If your physician prescribed oxygen for your sleep therapy, you can add the oxygen in one of two ways:

1. Obtain the optional oxygen adapter (part #7353D-601) from your home care provider and connect it to the outlet port on the SleepCube. Connect one end of the air supply tubing to the open end of the adapter and the other end of the air supply tubing to the mask.
2. If your mask is equipped with an oxygen port, your physician may request that you apply oxygen directly to your mask.

WARNING-OXYGEN USE

Oxygen supports combustion. To avoid possible physical injury, do not smoke while using this device with supplemental oxygen. Do not use this device near hot objects, volatile substances or sources of open flames.

Always turn on the SleepCube before turning on the oxygen source. Turn the oxygen source off before turning off the SleepCube.

This device is not suitable for use in the presence of a flammable anesthetic mixture with air or oxygen, or nitrous oxide.

At a fixed flow rate of supplemental oxygen, the inhaled oxygen concentration will vary depending on the pressure setting, patient breathing pattern, mask selection, and leak rate. This warning applies to most types of CPAP devices.

Never allow the oxygen source to run continuously while connected to the SleepCube if the device is not in use. If the SleepCube is not being used, turn off the oxygen flow.

CAUTION—*Oxygen is a prescription gas and should only be administered under the supervision of a physician.*

CAUTION—*The flow setting of the oxygen source must be specified by a physician.*

ACCESSORIES/REPLACEMENT ITEMS

Replacement Items

Air-inlet filter (4/pk)	DV51D-602
Optional fine particle filter (4/pk)	DV51D-603
Air supply tubing (6')	DV51D-629
Air supply plug	DV51D-604
Heater connection cover	DV51D-605
Oxygen adapter	7353D-601
Carrying case	DV51D-610

Power Items

AC power cords	
USA	DV51D-606
Europe, except U.K.	DV51D-607
U.K.	DV51D-608
Australia	DV51D-609
DC power cord	DV51D-619
DC battery clamp-on adapter	DV51D-696

NOTE- DC Battery clamp-on adapter requires DC power cord DV51D-619.

Masks Recommended For Use With The DeVilbiss SleepCube:

DeVilbiss FlexSet® Masks	9354 Series
DeVilbiss Serenity® Masks	9352 Series

Optional Heated Humidifier

DeVilbiss Heated Humidifier (cradle and chamber)	DV5HH
--	-------

SLEEPCUBE MAINTENANCE

DANGER

Electric shock hazard – Do not attempt to open or remove the enclosure; there are no user-serviceable components inside. If service is required, contact your home care provider for instructions on obtaining service of your device. Opening or attempting to service your device will void the warranty.

WARNING

If the operation of the device appears to deviate from normal operation or begins to produce unusual sounds, discontinue use and contact your home care provider for assistance.

If the device has been dropped, verify functionality of the device. If the device does not function or the enclosure is damaged, discontinue use and contact your provider.

If the device does not appear to be functioning properly refer to the Troubleshooting section for instructions.

Pressure Accuracy—The SleepCube is designed to provide accurate pressure performance from the factory without additional calibration. No routine service is required for a period of 2 years of use, provided the device is used in accordance with the manufacturer’s directions. After 2 years the device should be inspected.

NOTE—*Some countries require periodic servicing and calibration of this kind of medical device. Contact your home care provider for further information.*

Standard Air-Inlet Filter—Check the filter every 10 days and clean it if necessary. Replace the filter every 6 months or sooner if it is damaged. The standard filter is designed to filter particles larger than 3.0 microns in size. The SleepCube **MUST** have the standard filter in place during operation.

Optional Fine Particle Filter—Check the filter every 10 days, and replace it if it is dirty or damaged; otherwise replace the filter every 30 days. The optional fine particle filter is designed to filter particles as small as 0.3 microns in size and larger.

Filter Maintenance—Refer to Cleaning for instructions.

SLEEP-CUBE CLEANING

WARNING

To avoid electric shock, always unplug the power cord from the wall outlet power source before cleaning.

CAUTION—*Never rinse or place the device in water. Never allow liquids to get into or around any of the ports, switches or air filter; doing so will result in product damage. If this occurs, do not use the device. Contact your home care provider for service.*

Device Cabinet

1. Unplug the SleepCube from the power source. Wipe the enclosure with a clean, damp cloth every few days to keep the enclosure dust free.
2. Allow the device to dry completely before plugging it into a power source.

Filters

Air-Inlet Filter (Standard)

1. The air-inlet filter should be checked every 10 days and cleaned as necessary. The filter is located in the rear of the device.

CAUTION—*Proper filter function is important for the operation of the device and to protect the device from damage.*

2. To clean, remove the dark outer foam filter from the frame.
3. Wash the filter in a solution of warm water and dishwashing detergent, and rinse with water. Allow the filter to dry completely. If the filter (part #DV51D-602) is damaged, contact your home care provider for a replacement.
4. Ensure the filter is completely dry before installation and use.

Fine Particle Filter (Optional)

1. If you have an optional fine particle filter (part #DV51D-603), inspect this filter and replace it if it shows signs of dirt or damage. Do not wash.
2. Install the fine particle filter behind the standard filter.

NOTE—*If the fine particle filter is not installed first according to the directions, the life of the filter will be reduced requiring more frequent replacement.*

Tubing

The air supply tubing should be cleaned daily. Remove the tubing from the device and mask. Use a mild detergent and water to clean the inside of the tubing. Rinse the tubing and allow to air dry.

Mask and Headgear

Clean the mask and headgear according to the manufacturer’s instructions.

TROUBLESHOOTING

Electric Shock Hazard - Do not remove the device cabinet; there are no user-serviceable internal components. The cabinet should only be removed by a qualified DeVilbiss provider.

If you are having trouble with the operation of your DeVilbiss SleepCube, please refer to the table below for a list of problems with possible causes and remedies. If you are unable to resolve the problem by following the instructions, contact your home care provider for assistance.

Problem	Probable Cause	Remedy
Nothing shows on the display.	<ol style="list-style-type: none"> 1. Device is not plugged in or the power cord is not fully inserted. 2. There is no outlet power. 	<ol style="list-style-type: none"> 1a. Verify that power cord is firmly connected to the SleepCube and the power outlet. 1b. If you are using a DC power source make sure the cable connections are secure. Verify that your battery is charged. 2. Locate a suitable functioning power source.
The SleepCube does not start when breathing into the mask.	<ol style="list-style-type: none"> 1. The Auto-ON feature has been disabled. 2. There is no power to the device. 3. Breath is not deep enough for the Auto-ON feature to detect. 4. You are using a full face mask that has an anti-asphyxia valve. 5. The air supply port plug is missing or not fully inserted. 6. The humidifier chamber is not fully engaged into the cradle or is missing. 	<ol style="list-style-type: none"> 1. Use the ON/OFF button to start and stop the device. 2. Verify that the power cord is firmly connected to the device and the power outlet. 3. Breathe deeply in and out to start the SleepCube 4. Auto-ON may not work because your exhaled breath escapes through the open valve. Use the ON/OFF button to start and stop the device. 5. Make sure the air supply port plug is fully inserted into the unit. 6. Slide the humidifier chamber fully into the slot. Verify the latching lever snaps into place.

Problem	Probable Cause	Remedy
The airflow has stopped unexpectedly during use.	<ol style="list-style-type: none"> 1. The Auto-OFF feature has detected large airflow because of loosely fitting mask. 2. During use your mouth opens and you begin to mouth breathe. 	<ol style="list-style-type: none"> 1. Ensure you have a good mask fit; adjust the mask and headgear as needed. 2. Contact your home care provider for a chin strap or different mask which will help prevent mouth breathing. <p>If problem is not resolved with above remedies, disable Auto-OFF feature in Advanced Menu options.</p>
The display shows a device fault.	An error has occurred in the device and requires servicing.	Contact your home care provider for service.
Symptoms of OSA have recurred.	<ol style="list-style-type: none"> 1. Air filter may be occluded. 2. Sleep apnea condition has changed. 	<ol style="list-style-type: none"> 1. Clean or replace the air filter and relocate device away from drapes and other dusty surfaces. 2. Contact your physician.
Your skin becomes irritated where the mask touches your face.	<ol style="list-style-type: none"> 1. Headgear is too tight or improperly adjusted. 2. Your mask may not be properly sized or is not the most appropriate shape for you. 	<ol style="list-style-type: none"> 1. Loosen the headgear to reduce the contact pressure on your face. 2. Contact your physician or home care provider.
Dryness of throat or nose.	Inadequate humidity.	<ol style="list-style-type: none"> 1. Add a DeVilbiss DV5 Heated Humidifier. 2. Increase the heater setting on the humidifier.
Condensation of water is collecting in the hose causing a gurgling noise.	<ol style="list-style-type: none"> 1. Your humidification setting is too high. 2. The room temperature fluctuates from higher to lower levels through the night. 	Reduce the humidifier heater setting and/or increase the temperature of the room.
Air from flow generator seems to be too warm.	<ol style="list-style-type: none"> 1. The air filter(s) is/are dirty. 2. The air inlet port is blocked. 3. The room temperature is too high. 4. The SleepCube is located near a heater. 	<ol style="list-style-type: none"> 1. Refer to the Cleaning section for instructions on filter maintenance. 2. Unblock the air inlet. 3. Reduce the room temperature. 4. Locate the device away from the heat source so the air being drawn into the device is at room temperature. <p>If the problem persists contact your home care provider.</p>

SPECIFICATIONS

Problem	Probable Cause	Remedy
Nasal, sinus or ear pain, possibly with runny nose.	You may have a reaction to air flow pressure.	Discontinue use and contact your physician.

SPECIFICATIONS

Size	10.7 cm x 16.5 cm x 17.5 cm (4.2" H x 6.5" W x 6.9" D)
Weight	1.22 kg (2.7 lbs.)
Electrical Requirements AC	100-240V~, 50/60 Hz
Electrical Requirements DC	10.4 - 15 VDC, 5 Amps
Maximum Power Consumption	65 watts max from AC power source (flow generator only)
Pressure Range	3-20 cmH2O
Operating Temperature Range	5°C to 40°C (41°F to 104°F)
Operating Humidity Range	0% to 95% RH non-condensing
Operating Atmospheric Conditions	Sea level to 9,000 feet (2743 m)
Storage & Transportation Temperature Range	-40°C to +70°C (-40°F to 158°F)
Storage & Transportation Humidity Range	0% to 95% RH non-condensing
Maximum Limited Pressure	.20 cmH2O under normal use
Sound Pressure Level (tested per ISO 17510-1:2007)	26 dBA
Sound Power Level (tested per ISO 17510-1:2007)	38 dBA

Filter Specifications

- Standard Filter > 3.0 micron particles
- Optional Fine Particle Filter > 0.3 micron particles

Maximum Flowrate (per ISO 17510-1:2007)

	Test Pressures				
	3.0 cmH2O	7.0 cmH2O	12.0 cmH2O	16.0 cmH2O	20.0 cmH2O
Measured pressure at the patient connection port (cmH2O)	2.9	6.4	11.0	15.0	19.0
Average flow at the patient connection port (l/min)	83.8	126.8	137.4	130.0	113.8

Long Term Static Pressure Accuracy +/- 0.5 cmH2O

Dynamic Short-Term Pressure Accuracy with and without DV5HH Humidifier (per ISO 17510-1:2007)

Waveform	Volume [mL]	Breath Rate [min ⁻¹]	Set Pressure [cmH ₂ O]	Pressure Accuracy [Pk-Pk cmH ₂ O]
Sinusoidal cycle with I:E = 1:1	500	10, 15, 20	3.0	0.5
			7.0	0.5
			12.0	1.0
			16.0	1.0
			20.0	1.0

Equipment classification with respect to protection from electric shock Class II
 Degree of protection from electric shock Type BF Applied Part
 Degree of protection against ingress of liquids IPX1
 Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
 Mode of operation Continuous

Clinical Specifications for Digital Outputs

Mask Pressure 0 to 25.5 cmH₂O ± 0.5
 Patient Flow -127 to +127 L/min ± 10 L/min
 Leak Flow 0 to +127 L/min ± 10 L/min
 Estimated Tidal Volume 0 to 1023 mL

NOTE—Tidal Volume is an estimate. It is provided for trending purposes only.

DEVILBISS GUIDANCE AND MANUFACTURER'S DECLARATION

WARNING

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the Electromagnetic Compatibility [EMC] information provided in the accompanying documents.

Portable and Mobile RF Communications Equipment can affect Medical Electrical Equipment.

The equipment or system should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

NOTE– *The EMC tables and other guidelines provide information to the customer or user that is essential in determining the suitability of the Equipment or System for the Electromagnetic Environment of use, and in managing the Electromagnetic Environment of use to permit the Equipment or System to perform its intended use without disturbing other Equipment and Systems or non-medical electrical equipment.*

Guidance and Manufacturer's Declaration – Emissions All Equipment and Systems		
This device is intended for use in the electromagnetic environment specified below. The customer or user of this device should assure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Enforcement – Guidance
RF Emissions CISPR 11	Group 1	This device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B Radiated and Conducted Emissions	This device is suitable for use in all establishments including domestic, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonics IEC 61000-3-2	Class A	
Flicker IEC 61000-3-3	Complies	

Guidance and Manufacturer's Declaration – Emissions All Equipment and Systems			
This device is intended for use in the electromagnetic environment specified below. The customer or user of this device should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Floors should be wood, concrete or ceramic tile. If floors are synthetic, the relative humidity should be at least 30%.
Electrical Fast Transient/burst IEC 61000-4-4	±2kV on AC Mains	±2kV on AC Mains	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV Differential ±2kV Common	±1kV Differential ±2kV Common	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	>95% Dip for 0.5 Cycle 60% Dip for 5 Cycles 30% Dip for 25 Cycles >95% Dip for 5 Seconds	>95% Dip for 0.5 Cycle 60% Dip for 5 Cycles 30% Dip for 25 Cycles >95% Dip for 5 Seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or battery.

Guidance and Manufacturer's Declaration – Emissions All Equipment and Systems			
This device is intended for use in the electromagnetic environment specified below. The customer or user of this device should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Power Frequency 50/60Hz Magnetic Field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be that of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 Vrms from 150 kHz to 80 MHz	VI = 3 Vrms	Portable and mobile RF communications equipment should be separated from the device by no less than the recommended separation distances calculated/listed below: $D = (3.5/VI)\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	EI = 3V/m	$D = (3.5/EI)\sqrt{P}$ 80 to 800 MHz $D = (7/EI)\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum power rating in watts and D is the recommended separation distance in meters. Field strengths from fixed transmitters, as determined by an electromagnetic site survey, should be less than the compliance levels (VI and EI). Interference may occur in the vicinity of equipment containing a transmitter.
For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.			

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and this device. This device and system are NOT Life-Supporting

This device is intended for use in the electromagnetic environment in which radiated disturbances are controlled. The customer or user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF Communications Equipment and the device as recommended below, according to the maximum output power of the communications equipment.

Maximum Output Power (Watts)	Recommended Separation Distances for the device (meters)		
	150 kHz to 80 MHz	80 to 800MHz	800 MHz to 2.5 GHz
	$D=(1.1667)\sqrt{P}$	$D=(1.1667)\sqrt{P}$	$D=(2.3333)\sqrt{P}$
0.01	0.11667	0.11667	0.23333
0.1	0.36894	0.36894	0.73785
1	1.1667	1.1667	2.3333
10	3.6894	3.6894	7.3785
100	11.667	11.667	23.333

For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

DECLARATION OF CONFORMITY

Manufacturer: DeVilbiss Healthcare LLC
Address: 100 DeVilbiss Drive
Somerset, Pennsylvania 15501-2125 USA

Product Designation: CPAP UMDNS 11-001

Type, Model: DeVilbiss IntelliPAP or SleepCube DV5x Series CPAP
– Models DV51 Standard CPAP & DV54 AutoAdjust
DeVilbiss Heated Humidifier
– Model DV5HH

We herewith declare that the above mentioned product complies with the requirements of EC Directive 93/42/EEC + A2007 Annex II, 3

Class: IIa, Rule 9

Quality System
Standards Applied: ISO13485:2003

Notified Body: TÜV NORD Cert GmbH

Safety Standards Applied: IEC 60601-1:1988 + AI 1995
IEC 60601-1-4:2000
ISO 14971:2007
EN ISO17510-1:2007
IEC 60601-1- 6:2004

EMC Compliance To: IEC 60601-1-2:2001 + AI 2004

Authorized Representative: DeVilbiss Healthcare GmbH
Kamenzer Straße 3
68309 Mannheim
Germany
+49 621-178-98-0

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DeVilbiss Healthcare Pty. Limited
DeVilbiss Healthcare SAS
DeVilbiss Healthcare LLC • 100 DeVilbiss Drive • Somerset, PA 15501 • USA